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Name and Address of Manufacturer and EU Authorised Representative

	Micronclean Ltd Roman Bank, Skegness, Lincs, PE25 1SQ, UK Tel: +44 (0)1754 767377 www.micronclean.com
	Advena Ltd. Tower Business Centre, 2nd Flr. Tower Street, Swatar. BKR 4013 Malta

Sterility

The devices are supplied in a sterile state and the method of sterilisation is gamma irradiation. The products should not be used if the primary packaging is open or damaged as the sterility of the devices could be compromised.

UKCA / CE Marking

The devices are CA marked as per the Medical Devices Regulations 2002 (UK MDR 2002) and CE marked as per the Medical Device Regulation (EU) 2017/745. Each mark is accompanied by the number of the applicable approved / notified body. The marks appear as below:

**Intended Purpose**

The devices are single use only. Risks of re-use of the devices include:

- Syringes and caps would have to be returned to a compounding unit for re-filling in a non-sterile state with a risk of microbial contamination
- Incompatibility of one drug product with another leading to reduced drug efficacy
- The syringes are not qualified for any decontamination procedures; decontamination to allow re-use could introduce new risks not known to, or quantifiable by, Micronclean

The intended use of the syringes and caps is the manual administration of pharmaceutical products to patients (the syringe having been pre-filled with the pharmaceutical product in a compounding pharmacy). Micronclean do not supply needles for injection, these are supplied by a healthcare facility at the point of use. The syringes are not intended for use in automated syringe pumps.

The medical conditions which can be treated using the syringes and caps supplied by Micronclean are determined by healthcare professionals. Micronclean do not make any specific clinical claims for the devices.

The syringes and caps have luer lock fittings and can be used in conjunction with other luer lock components and devices, as determined by healthcare professionals.

Syringe caps are supplied for use to cap a pre-filled syringe. The cap twists onto the syringe tip and protects the tip prior to use, prevents leakage from the tip and forms a sterile barrier with the syringe to prevent contamination of the syringe and its contents.

The syringes should not be used for intra-ocular purposes. This is due to the risk of silicone oil depositing in the vitreous layer of the eye.

The devices should not be used if damaged or contaminated.

The syringe manufacturer BD does not support the use of their syringes as drug storage containers and recommend that the syringes are used for general drawing up and aspiration only. Micronclean recognise that many end users use the syringes as a drug storage container and support this use provided that the user has carried out a suitable qualification (or has access to qualification data) eg drug stability studies, and meets all applicable GMP and regulatory requirements when compounding the drug.

Sealed pouches have a lifetime of three years. Pouches should be opened in a clean area using suitable aseptic techniques. Devices should be used as soon as possible following the opening of the pouch.

Performance Characteristics

The performance characteristics of the devices are:

- Products are sterile (sterility assurance level 10-6)
- Sealed pouches have a three year shelf-life
- Products are non-invasive
- Under qualified conditions, the outer surface of the packaging is suitable for sanitisation by VHP and alcohol methods. Other sanitisation methods may be appropriate if qualified by the end user
- Syringes deliver an accurate volume of liquid (tolerance as per ISO 7886-1)

Storage Conditions

Store in a dry, well ventilated area and protect from direct sunlight and heat sources.

Device Disposal

If required by local protocols and procedures, syringes and syringe caps should be decontaminated prior to disposal to remove any drug residues.

Syringe barrels and plungers and syringe caps are generally suitable for recycling, subject to the recycling contractor's requirements. The syringe stopper (aka bung or gasket) is generally not suitable for recycling.

Decontaminated syringes and syringe caps are suitable for landfill. Where non-CA marked drawing up needles are used in conjunction with the syringes, these should be disposed of in a suitable sharps container.

Serious Incidents

In the event of a serious incident (ie one which led to, or could have led to the death or serious illness of a patient) the user should inform Micronclean without delay, on **+44 1754 767377**. The national competent authority in your country should also be informed.