



## PRODUCT SPECIFICATION

### POLYESTER POUCH WIPE STERILE AND NON-STERILE

VeriGuard **3**

#### PRODUCT INFORMATION

##### PRODUCT DESCRIPTION

VeriGuard 3 is a pre-saturated Polyester wipe designed for use in high grade Cleanrooms as a hard surface wipe.

##### PRODUCT FEATURES

100% continuous filament pouch wipe pre-saturated with IPA 70/ 30 WFI supplied in a resealable pouch.

Each wipe is 300mm x 300mm supplied either non-folded or V-Folded to suit your requirements.

##### MANUFACTURING

Primary packaged in polythene bags in ISO 14644 Class 5 (at rest) cleanroom.

##### DISINFECTANT SPECIFICATION

Contains: 70% Propan-2-ol (#67-63-0) and 30% Water for Injection grade water.

##### PACKAGING

Packed in cleanroom compatible pouch containing either 50 or 15 wipes.  
Double Bagged (sterile only).

##### QUALITY ASSURANCE

Product supplied under BS EN ISO 9001.

##### SHELF LIFE

2 year from date of manufacture.

##### STORAGE CONDITIONS

Store in a cool, dry, well ventilated area and protect from direct sunlight heat and naked flames.

## MATERIALS & PERFORMANCE

Performance Data			
Product Code	Weight	Size	Thickness
MCZNLA50N-A/I	130±10g/m <sup>2</sup>	300mm x 300mm	
MCZSLA50S-A/I	130±10g/m <sup>2</sup>	300mm x 300mm	
WNVG06030-A/I	120 + 5 gsm	300mm x 300mm	0.52±0.04mm
WSVG06030-A/I	120 + 5 gsm	300mm x 300mm	0.52±0.04mm

## PRODUCT INFORMATION

### DISINFECTANT PERFORMANCE

SUSPENSION TESTS					
Protocol	Organism	Contact Time (min)	Conditions	Log reduction to pass	Log reduction
BS EN 1040	S. aureus	5 min	Dirty	Log 5.0	> 5.3
	P. aeruginosa	5 min	Dirty	Log 5.0	> 5.4
BS EN 1275	C. albicans	5 min	Dirty	Log 4.0	> 4.4
BS EN 1276	S. aureus	5 min	Clean	Log 5.0	> 5.3
	E. coli	5 min	Clean	Log 5.0	> 5.3
	P. aeruginosa	5 min	Clean	Log 5.0	> 5.4
	E. hirae	5 min	Clean	Log 5.0	> 5.3
BS EN 1650	C. albicans	5 min	Clean	Log 4.0	> 4.4

  

SURFACE TESTS					
Protocol	Organism	Contact Time (min)	Conditions	Log reduction to pass	Log reduction
BS EN 13697	S. aureus	5 min	Clean	Log 4.0	6.0
	E. coli	5 min	Clean	Log 4.0	> 4.8
	P. aeruginosa	5 min	Clean	Log 4.0	> 4.1
	E. hirae	5 min	Clean	Log 4.0	> 5.5
	C. albicans	5 min	Clean	Log 3.0	> 4.0
	A. niger	5 min	Clean	Log 3.0	3.7

Please note: Data shown above are just typical values.

## ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
MCZNLA50N		IPA Laundered Polyester Non-Folded Pouch Wipe Non-Sterile	300 x 300mm	50	5 x 50	390 x 390 x 190mm	8.0kg
	MCZSLA50S	IPA Laundered Polyester Non-Folded Pouch Wipe Sterile	300 x 300mm	50	5 x 50	390 x 390 x 190mm	8.0kg
WNVG06030		IPA Polyester V-Folded Pouch Wipe Non-Sterile	300 x 300mm	15	30 x 15	400 x 400 x 395mm	15.78kg
	WSVG06030	IPA Polyester V-Folded Pouch Wipe Sterile	300 x 300mm	15	30 x 15	400 x 400 x 395mm	15.78kg

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please contact us.