

CONSUMABLES

INTERNATIONAL BROCHURE :: 2023





CONTENTS

INTRODUCTION	01	GOGGLES	29
REGIONAL ORDERING AND LABELS	02	INDIVIDUALLY PACKED GOGGLE NON-STERILE	30
ALCOHOLS	03	CLOTHING	31
TRIGGER SPRAY FEATURES	04	DISPOSABLE COVERALL WITH INTEGRAL HOOD	32
TRIGGER SPRAY	05	DISPOSABLE COVERALL WITH INTEGRAL HOOD AND FEET STERILE	33
RFU 5L	06	DISPOSABLE LABCOAT NON-STERILE	34
IPA POLYPROPYLENE TUB WIPE	07	CLEANROOM POLYESTER SOCK	35
IPA AND IMS POLYCELLULOSE POUCH WIPE	08	DISPOSABLE FACEMASKS, BEARD SNOOD AND MOB CAP NON-STERILE	36
IPA POLYESTER POUCH WIPE	09	MATS	37
IMS SILVER PACK POLYETHYLENE/POLYPROPYLENE PRE-SATURATED WIPE	10	ACCESS MAT MEDIUM TACK NON-STERILE	38
NEUTRAL DETERGENT	11	ACCESS MAT LOW TACK NON-STERILE	39
DELTA TRIGGER SPRAY	12	RECYCLING INFORMATION	40
DELTA AND DI 5L RFU	13		
DELTA POLYCELLULOSE POUCH WIPE	14		
DISPOSABLE MOP SYSTEM	15		
ISOLATOR BONNET AND TOOL	16		
IMS DISPOSABLE MICROFIBRE MOP STERILE	17		
DELTA DISPOSABLE MICROFIBRE MOP STERILE	18		
DRY DISPOSABLE MICROFIBRE MOP STERILE	19		
WASTE BAGS, MOP HEADS, TROLLEY AND OTHER ACCESSORIES	20		
DRY WIPES	21		
POLYCELLULOSE WIPE STERILE	22		
POLYCELLULOSE WIPE NON-STERILE	23		
POLYESTER WIPE STERILE	24		
POLYESTER WIPE NON-STERILE	25		
MICRODENIER WIPE NON-STERILE	26		
POLYCELLULOSE ROLL WIPE NON-STERILE	27		
LAUNDERED DUAL LAYER POLYESTER WIPE STERILE	28		



INTRODUCTION

"At Micronclean, our passion is to be the first to develop new technological solutions that change the shape of the markets we serve creating quality and efficiency for our customers."

For over 40 years Micronclean has focused on contamination control products for companies in the pharma, healthcare, microelectronics, aerospace and semiconductor industries. As a result, we are incredibly proud of our cleanroom consumables range which we have sold to over 30 countries around the world in the last 5 years.

The 2023 edition of our Consumables Catalogue has been updated to include new products such as the Autoclavable Goggles, alongside recycling information for all our products.

A large percentage of our range is manufactured in the UK under cleanroom conditions. For those products we do not manufacture we have full control of the supply chain and set the specification. Having the capability to manufacture means we are able to offer bespoke products and maintain sufficient stock levels to meet customer demand and our customers are able to audit our facilities.

For more information on any of our product lines please visit our website www.micronclean.com or contact us directly on cs@micronclean.co.uk or +44 (0) 1754 767377

Michelle Simpson
Product Director

REGIONAL ORDERING AND LABELS

MICRONCLEAN UPDATES ITS PRODUCT LABELS TO DEAL WITH DEMAND FROM EXPORT MARKETS.

Following increasing demand for our market leading disinfectant and cleaning products, we have had to enhance our existing product labels to ensure that we remain compliant.

As a result, we have generated an enhanced label which enables us to remain compliant with labeling regulations and allows us to carry up to three languages in addition to English.

So, the result of these changes is that the product code for our existing CleanGuard 1 IPA 500ml Trigger Spray Sterile will migrate from the current ZSAL64499 to ZSAL64499-A for region A orders, ZSAL64499-B for region B orders and so on. Note this regionalisation only affects products with active chemicals.

The new look labels with the associated regional codes for products will be rolled out progressively across our whole consumable range.

Regional Labels in Use

For a UK consumable customer, the default option will always be that they order region A products and they will normally receive region A products. However, as the English language is present on all the labels and labelling regulations are common across region A-E, we might on occasion substitute the same product with a region B label for instance.

For a Micronclean customer in Italy for instance, the product should always be ordered for region B and to comply with the local language regulations we will supply the product with region B labelling.

Although English is present on all the labels, for regions G and I the products are destined for Australia and India respectively. This is because they have to conform to a different set of labeling regulations in each of these countries and as such these products can only be sold and used in these countries.

Please note if you have any questions please do not hesitate to contact us.


REGIONS	A	UK Ireland Germany Austria
	B	UK Spain Portugal Italy
	C	UK France Belgium Holland
	D	UK Sweden Denmark Norway
	E	UK Greece Poland
	G	UK Australia
	I	UK India

ALCOHOLS



ALCOHOL TRIGGER SPRAY FEATURES

- Long trigger handle for ease of use
- Mechanically sealed trigger head
- Contains 70% alcohol blended with WFI quality water
- No bag in the bottle for easy viewing of liquid level through transparent bottle
- Contents filtered to 0.2µm at point of fill
- Ergonomic low centre of gravity bottle design





RECYCLING INFORMATION



+ info see page 40

ALCOHOL TRIGGER SPRAY

CLEANGUARD 1 & 2

USAGE

Designed for use in Cleanrooms
Cleanroom transfer process
Disinfection of hard surfaces
Use within Isolators / LAFs

BENEFITS

No bag in the bottle
70% alcohol blended with WFI quality water
Filtered to 0.2µm at point of fill
Mechanically sealed trigger head
Low centre of gravity bottle design
Transparent bottle
Long trigger handle, ease of use

FEATURES

Trigger spray with spray / off feature
2 year shelf life
Gamma irradiation stable
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
ZSAL64499	IPA 500ml Trigger Spray Sterile	500ml	Each	24 x 500ml
ZSAL64950	IPA 950ml Trigger Spray Sterile	950ml	Each	15 x 950ml
ZNAL64950	IPA 950ml Trigger Spray Non-Sterile	950ml	Each	15 x 950ml
ZSAL65499	IMS 500ml Trigger Spray Sterile	500ml	Each	24 x 500ml
ZSAL65950	IMS 950ml Trigger Spray Sterile	950ml	Each	15 x 950ml
ZNAL65950	IMS 950ml Trigger Spray Non-Sterile	950ml	Each	15 x 950ml

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



RECYCLING INFORMATION

+ info see page 40



IPA AND IMS 5 Litre RFU CLEANGUARD 1 & 2

USAGE

Disinfection of hard surfaces

BENEFITS

70% Alcohol blended with WFI quality water
Bulk container for larger surface areas
Easy to use with bucket mopping systems
Filtered to 0.2µm at point of fill

FEATURES

Tamper evident PP vented cap
Gamma irradiation stable
2 year shelf life
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
ZSAL64005	IPA 5 Litre RFU Sterile	5L	Each	2 x 5L
ZSAL65005	IMS 5 Litre RFU Sterile	5L	Each	2 x 5L

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



RECYCLING INFORMATION



+ info see page 40

IPA POLYPROPYLENE TUB WIPE

VERIGUARD 1

USAGE

Hard surface wiping particularly suited to outer areas

BENEFITS

Designed for low grade cleanroom areas
Designed for high volume usage
Long lasting alcohol retention
Re-fillable pouch available

FEATURES

Continuous filament thermobonded nonwoven polypropylene wipe
18gsm wipe weight
2 years shelf life
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
WSVG05020	IPA Polypropylene Tub Wipe Sterile	200mm x 200mm	200	12 x 200
WSVG01020	IPA Polypropylene Tub Wipe Refill Sterile	200mm x 200mm	200	20 x 200
WNVG05020	IPA Polypropylene Tub Wipe Non-Sterile	200mm x 200mm	200	12 x 200

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



RECYCLING INFORMATION

+ info see page 40



IPA AND IMS POLYCELLULOSE POUCH WIPE VERIGUARD 2

USAGE

Hard surface wiping particularly suited to transfer processes

BENEFITS

Increased alcohol retention from re-sealable pouch
High absorbency
Strong multifunctional wipe
Small pack size available for small batch manufacturing, large pack size suitable for large batch manufacturing

FEATURES

White hydro-entangled Polycellulose wipe
68gsm wipe weight
2 years shelf life
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
WSVG04023	IPA Polycellulose C-folded Pouch Wipes Sterile	230mm x 230mm	100	8 x 100
WSCW92009	IPA Sessional Polycellulose C-folded Pouch Wipes Sterile	230mm x 230mm	15	40 x 15
WNVG04023	IPA Polycellulose C-folded Pouch Wipes Non-Sterile	230mm x 230mm	100	8 x 100
WSVG09023	IMS Polycellulose C-folded Pouch Wipes Sterile	230mm x 230mm	100	8 x 100
WSCW91009	IMS Sessional Polycellulose C-folded Pouch Wipes Sterile	230mm x 230mm	15	40 x 15

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

RECYCLING INFORMATION



+ info see page 40



IPA POLYESTER POUCH WIPE

VERIGUARD 3

USAGE

Cleanroom hard surface wipe
Use in Isolators / Laminar flow cabinets

BENEFITS

Increased alcohol retention from re-sealable pouch
Laser cut edges
Strong multifunctional wipe
Small pack size available for small batch manufacturing, medium pack size suitable for larger batch manufacturing

FEATURES

100% continuous filament polyester wipe
2 years shelf life
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
WSVG06030	IPA Sessional Polyester V-folded Pouch Wipe Sterile	300mm x 300mm	15	30 x 15
MCZSLA50S	IPA Polyester Laundered Non-folded Pouch Wipes Sterile	300mm x 300mm	50	5 x 50
WNVG06030	IPA Sessional Polyester V-folded Pouch Wipe Non-Sterile	300mm x 300mm	15	30 x 15
MCZNLA50N	IPA Polyester Laundered Non-folded Pouch Wipes Non-Sterile	300mm x 300mm	50	5 x 50

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



RECYCLING INFORMATION

+ info see page 40



IMS SILVER PACK POLYETHYLENE / POLYPROPYLENE PRE-SATURATED WIPE VERIGUARD 9

USAGE

Cleanroom hard surface wipe

BENEFITS

- Increased alcohol retention from re-sealable pouch
- Designed for large batch manufacturing
- Blended with de-ionised water (DI)
- Long lasting alcohol retention

FEATURES

- Nonwoven interleaved wipes
- 2 years shelf life
- Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
WSCW97ZZZ	IMS Silver Pack Polyethylene / Polypropylene Z-folded Wipes Sterile	200mm x 220mm	200	8 x 200

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

NEUTRAL DETERGENT





DELTA TRIGGER SPRAY

CLEANGUARD 5

USAGE

Cleaning of hard surfaces

Annex 1 states that for disinfection to be effective it must be a two-stage process. The first stage should be a process to remove surface contamination and the second stage should be disinfection using a rotational system based on two types of disinfectants which have different modes of action

Micronclean has developed Delta as a sterile neutral detergent to remove surface contamination

BENEFITS

No bag in the bottle
Blended with WFI quality water
Filtered to 0.2µm at point of fill
Mechanically sealed trigger head
Low centre of gravity bottle design
Transparent bottle
Long trigger handle, ease of use

FEATURES

Trigger spray with spray / off feature
2 year shelf life
Gamma irradiation stable
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
ZSND05950	Delta 950ml Trigger Spray Sterile	950ml	Each	15 x 950ml

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

RECYCLING INFORMATION



+ info see page 40



DELTA AND DI 5L RFU CLEANGUARD 5 AND 6

USAGE

Cleaning of hard surfaces and floors

Annex 1 states that for disinfection to be effective it must be a two-stage process. The first stage should be a process to remove surface contamination and the second stage should be disinfection using a rotational system based on two types of disinfectants which have different modes of action

Micronclean has developed Delta as a sterile neutral detergent to remove surface contamination.

BENEFITS

Blended with WFI quality water

Bulk container for larger surface areas

Easy to use with bucket mopping systems

Filtered to 0.2µm at point of fill

FEATURES

Tamper evident PP vented cap

Gamma irradiation stable

2 year shelf life

Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
ZSND05005	Delta 5 Litre RFU Sterile	5L	Each	2 x 5L
ZSDI82005	DI Water 5 Litre RFU Sterile	5L	Each	2 x 5L
ZNDI82005	DI Water 5 Litre RFU Non-sterile	5L	Each	2 x 5L

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



RECYCLING INFORMATION

+ info see page 40



DELTA POLYCELLULOSE POUCH WIPE

VERIGUARD 2

USAGE

Hard surface wiping

Annex 1 states that for disinfection to be effective it must be a two-stage process. The first stage should be a process to remove surface contamination and the second stage should be disinfection using a rotational system based on two types of disinfectants which have different modes of action

Micronclean has developed Delta as a sterile neutral detergent to remove surface contamination

BENEFITS

- Increased liquid retention from re-sealable pouch
- Designed for small batch manufacturing
- High absorbency


FEATURES

- White hydro-entangled Polycellulose wipe
- 68gsm wipe weight
- 2 years shelf life
- Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
WSN05009	Delta Sessional Polycellulose C-folded Pouch Wipe Sterile	230mm x 230mm	15	10 x 15

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

DISPOSABLE MOP SYSTEM





ISOLATOR BONNET AND TOOLS

PUREGUARD 6 AND 7

USAGE

Designed specifically to help clean the inside of isolators, RABs and LAF cabinets


BENEFITS

Universal fit
Dual layer polyester guarantees high levels of absorbency
Lightweight isolator tool available in 2 handle lengths with 360° rotation

FEATURES

Low linting, dual-layer 100% continuous filament polyester for increased absorbency
250gsm bonnet weight
Supplied sterile one per bag ready for use
2 year sterile shelf life
Handle supplied Non-Sterile ready for sterilisation

Product Code	Description	Size	Pack Qty	Case Qty
ISED0140	Dry Dual Layer Polyester Laundered Isolator Bonnet Sterile	One Size	EA	40
ZNSS01035	Stainless Steel Isolator Cleaning Tool 35cm Non-Sterile	One size	EA	1
ZNSS01061	Stainless Steel Isolator Cleaning Tool 61cm Non-Sterile	One size	EA	1



IMS DISPOSABLE MICROFIBRE MOP STERILE PUREGUARD 2

USAGE

Designed for use on cleanroom floors

BENEFITS

Uniformly dosed to achieve a coverage of 15m²
No need for in cleanroom mixing or dosing
2 pocket design enables a hands-free and secure connection to the mop frame
Eliminates the risk associated with exposed chemicals and soiled water

FEATURES

Proven 2 year shelf life
Individually packed mops are suitable for smaller cleanrooms
Multi packed mops are suitable for larger cleanrooms
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
MSME24016	IMS Individually Packed Disposable Microfibre Mop Sterile	440mm x 140mm	1	30 Packs of 1
MSME21016	IMS Multi Packed Disposable Microfibre Mop Sterile	440mm x 140mm	10	3 Packs of 10

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



RECYCLING INFORMATION

+ info see page 40



DELTA DISPOSABLE MICROFIBRE MOP STERILE

PUREGUARD 5

USAGE

Annex 1 states that for disinfection to be effective it must be a two-stage process. The first stage should be a process to remove surface contamination and the second stage should be disinfection using a rotational system based on two types of disinfectants which have different modes of action

Micronclean has developed Delta as a sterile neutral detergent to remove surface contamination.

BENEFITS

Uniformly dosed to achieve a coverage of 15m²

No need for in cleanroom mixing or dosing

2 pocket design enables a hands-free and secure connection to the mop frame

Eliminates the risk associated with exposed chemicals and soiled water

FEATURES

Proven 2 year shelf life


Individually packed mops are suitable for smaller cleanrooms

Multi packed mops are suitable for larger cleanrooms

Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
MSMN24016	Delta Individually Packed Disposable Microfibre Mop Sterile	440mm x 140mm	1	30 Packs of 1
MSMN21016	Delta Multi Packed Disposable Microfibre Mop Sterile	440mm x 140mm	10	3 Packs of 10

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



DRY DISPOSABLE MICROFIBRE MOP STERILE

PUREGUARD 6

USAGE

Designed for use on cleanroom floors

BENEFITS

2 pocket design enables a hands-free and secure connection to the mop frame

FEATURES

Proven 2 year shelf life
Individually packed mops are suitable for smaller cleanrooms
Multi packed mops are suitable for larger cleanrooms
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
MSMD24016	Dry Individually Packed Disposable Microfibre Mop Sterile	One size	1	30 Packs of 1
MSMD21016	Dry Multi Packed Disposable Microfibre Mop Sterile	One Size	10	4 Packs of 10

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

Please see individual product specifications



RECYCLING INFORMATION



WASTE BAGS, MOP HEADS, TROLLEY AND OTHER ACCESSORIES

PUREGUARD 7 AND 8


BENEFITS

Designed to work alongside the PureGuard Mopping System


Product Code	Description	Size	Pack Qty	Case Qty
BSP1C1012	Polyethylene Bag Sterile	One Size	10	12 Packs of 10
MTSS01001	Stainless Steel Standard Mop Trolley Non-Sterile	One size	Each	1
MFSS01001	Stainless Steel Lightweight Mop Frame 40cm Non-Sterile	One Size	Each	1
MFSS79226	Stainless Steel Shrouded Locking Mop Frame 40cm Non-Sterile	One Size	Each	1
MHSS02001	Stainless Steel Fixed Mop Handle 142cm Non-Sterile	One Size	Each	1
MHSS02002	Aluminium Adjustable Mop Handle 100cm-180cm Non-Sterile	One Size	Each	1
MTSS01001	Stainless Steel Standard Mop Trolley Non-Sterile	One Size	Each	1

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the tables on page.

DRY WIPES







RECYCLING INFORMATION

+ info see page 41



POLYCELLULOSE WIPE STERILE

WHITEGUARD 2

USAGE

Hard surface wipe
Transfer of product into cleanrooms

BENEFITS

Excellent absorbency

FEATURES

Low linting, multipurpose wipes
68gsm wipe weight
Hydroentangled 55% Polyester, 45% Cellulose blend

Product Code	Description	Size	Pack Qty	Case Qty
WSWG04009	Polycellulose Non-Folded Wipe Sterile	230mm x 230mm (9" x 9")	100	10 Packs of 100
WSWG06009	Polycellulose Non-Folded Wipe Sterile	230mm x 230mm (9" x 9")	300	8 Packs of 300
WSWG04012	Polycellulose Non-Folded Wipe Sterile	300mm x 300mm (12" x 12")	100	10 Packs of 100
WSCS93009	Polycellulose Non-Folded Wipe Sterile	230mm x 230mm (9" x 9")	300	8 Packs of 300
WSCS93012	Polycellulose Non-Folded Wipe Sterile	300mm x 300mm (12" x 12")	10 x 10	17 Packs of 100



+ info see page 41



POLYCELLULOSE WIPE NON-STERILE WHITEGUARD 2

USAGE

Hard surface wipe
Transfer of product into cleanrooms




BENEFITS




Excellent absorbency

FEATURES

Low linting, multipurpose wipes
68gsm wipe weight
Hydroentangled 55% Polyester, 45% Cellulose blend

Product Code	Description	Size	Pack Qty	Case Qty
WNWG04009	Polycellulose Non-Folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	100	10 Packs of 100
WNWG06009	Polycellulose Non-Folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	300	16 Packs of 300
WNWG04012	Polycellulose Non-Folded Wipe Non-Sterile	300mm x 300mm (12" x 12")	100	20 Packs of 100
WNCS93009	Polycellulose Non-Folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	300	8 Packs of 300
WNCS93012	Polycellulose Non-Folded Wipe Non-Sterile	300mm x 300mm (12" x 12")	100	17 Packs of 100





RECYCLING INFORMATION

+ info see page 41



POLYESTER WIPE STERILE

WHITEGUARD 3

USAGE

Hard surface wipe
Transfer of product into cleanrooms

BENEFITS

Ultra low linting wipe

FEATURES

Laser sealed edges to prevent linting
100% knit continuous filament polyester
Medium levels of absorbency

Product Code	Description	Size	Pack Qty	Case Qty
WSTG05009	Polyester Non-Folded Wipe Sterile	230mm x 230mm (9" x 9")	150	10 packs of 150
WSTG03012	Polyester Non-Folded Wipe Sterile	300mm x 300mm (12" x 12")	50	12 packs of 50
WSPW80009	Polyester Laundered Non-Folded Wipe Sterile	230mm x 230mm (9"x9")	150	10 packs of 150
MCDELN50S	Polyester Laundered Single Layer Non-Folded Wipe Sterile	300mm x 300mm (12"x12")	50	10 packs of 50

RECYCLING INFORMATION



+ info see page 41



POLYESTER WIPE NON-STERILE WHITEGUARD 3

USAGE

Hard surface wipe
Transfer of product into cleanrooms



BENEFITS



Ultra low linting wipe



FEATURES

Laser sealed edges to prevent linting
100% knit continuous filament polyester
Medium levels of absorbency

Product Code	Description	Size	Pack Qty	Case Qty
WNTG04009	Polyester Non-Folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	100	10 packs of 100
WNTG05009	Polyester Non-Folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	150	10 packs of 150
WNTG05012	Polyester Non-Folded Wipe Non-Sterile	300mm x 300mm (12" x 12")	150	10 packs of 150
WNTG03012	Polyester Non-Folded Wipe Non-Sterile	300mm x 300mm (12" x 12")	50	10 packs of 50







RECYCLING INFORMATION

+ info see page 41



MICRODENIER WIPE NON-STERILE

WHITEGUARD 5

USAGE

Hard surface wipe
Excellent wipe for lens and screen cleaning

BENEFITS

Ultra low linting wipe

FEATURES

100% continuous filament nylon/polyester
wipe

Product Code	Description	Size	Pack Qty	Case Qty
WNMD04009	Microdenier Non-folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	100	10 Packs of 100



RECYCLING INFORMATION



+ info see page 41

POLYCELLULOSE ROLL WIPE NON-STERILE

WHITEGUARD 7

USAGE

For use in support areas
Hard surface wipe




BENEFITS




Excellent absorbency

FEATURES

Low linting, multipurpose wipes
Log roll of 400 wipes

Product Code	Description	Size	Pack Qty	Case Qty
WNPW81010	Polycellulose Roll Wipe Non-Sterile White	254mm x 381mm	400 sheets/roll	2 x 1 Roll
WNPB81010	Polycellulose Roll Wipe Non-Sterile Blue	254mm x 381mm	400 sheets/roll	2 x 1 Roll





RECYCLING INFORMATION

+ info see page 41



LAUNDERED DUAL LAYER POLYESTER WIPE STERILE WHITEGUARD 11

USAGE

Hard surface wipe
Transfer of product into cleanrooms

BENEFITS


Ultra low linting wipe

FEATURES

Dual layer 100% continuous filament polyester
Laser sealed edges to prevent linting
Medium levels of absorbency

Product Code	Description	Size	Pack Qty	Case Qty
MHPPELN50S	Dual Layer Polyester Laundered Non-folded Wipe Sterile	300mm x 300mm (12" x 12")	50	4 Packs of 50

GOGGLES





RECYCLING INFORMATION

+ info see page 41



INDIVIDUALLY PACKED GOGGLE NON-STERILE VISIONGUARD 1

USAGE

Specifically designed for use in
Cleanroom Environments

BENEFITS


Fully compatible with the requirements of
Annex 1


FEATURES

Supplied in a Non-Sterile condition but
packed in an Autoclave bag ready for
immediate sterilisation
40 process lifespan for goggles without
significant fogging
Lightweight, supplied with adjustable head
strap and easy-to-use snap-on clip
Indirect ventilation to top of goggle to prevent
fogging and protect from splashes
Compatible with most prescription glasses
EN 166 3 4 BT Compliant
CE and UKCA Marked

Product Code	Description	Size	Pack Qty	Case Qty
ONACL0112	Individually Packed Goggle Non-Sterile	One Size	1	12 Packs of 1

CLOTHING





+ info see page 41



DISPOSABLE COVERALL WITH INTEGRAL HOOD

SUREGUARD 1

USAGE

Designed for cleanroom use

BENEFITS

CE Marked to PPE Category III Type 5 & 6
Anti-static

FEATURES

Single use cleanroom coverall
CE and UKCA Marked.
Polypropylene with polyethylene laminate fabric
Integrated hood
Latex free elastic at hood, wrists and ankles
Available in Sterile and Non-Sterile Forms

Product Code	Description	Size	Pack Qty	Case Qty
DSWH04	Disposable Coverall with Integral Hood Sterile	SML - 2XL	1	25
DNWH04	Disposable Coverall with Integral Hood Non-Sterile	SML - 2XL	1	25



RECYCLING INFORMATION



+ info see page 41

DISPOSABLE COVERALL WITH INTEGRAL HOOD AND FEET STERILE

SUREGUARD 3

USAGE

Designed for cleanroom use


BENEFITS

CE Marked to PPE Category III Type 5 & 6
Anti-static

FEATURES

Single use cleanroom coverall
CE and UKCA Marked.
Polypropylene with polyethylene laminate fabric
Integrated hood
Latex free elastic at hood, wrists and ankles

Product Code	Description	Size	Pack Qty	Case Qty
DSWH21	Disposable Coverall with Integral Hood and Feet Sterile	SML - 2XL	1	25



+ info see page 41



DISPOSABLE LABCOAT NON-STERILE SUREGUARD 2

USAGE

Designed for laboratory and cleanroom use

BENEFITS

CE Marked to PPE Category III Type 5 & 6
Anti-static

FEATURES

Single use cleanroom labcoat
CE and UKCA Marked
Polypropylene with polyethylene laminate fabric
Latex free elastic at wrist

Product Code	Description	Size	Pack Qty	Case Qty
DNWH10	Disposable Labcoat Non-Sterile	SML - 2XL	1	25



RECYCLING INFORMATION



+ info see page 41

CLEANROOM POLYESTER SOCK

TUFFGUARD 1

FEATURES

Specifically developed for use in Cleanrooms to help with compliance with Annex 1

BENEFITS

Sock shaped for wearer comfort

FEATURES

White sock shaped
100% Polyester Socks
One Size

Product Code	Description	Size	Pack Qty	Case Qty
SNWH31ZZZ	Polyester Cleanroom Sock	One Size	50 (Single)	50



RECYCLING INFORMATION

+ info see page 41



FACEGUARD 3
FACEMASK WITH HEADLOOPS



FACEGUARD 1
FACEMASK WITH TIES



FACEGUARD 5 AND 6
BEARD SNOOD AND MOP CAP



FACEGUARD 2
FACEMASK WITH EARLOOPS

DISPOSABLE FACEMASK FACEGUARD 1, 2 AND 3 BEARD SNOOD AND MOB CAP NON-STERILE FACEGUARD 5 AND 6

FEATURES

DISPOSABLE FACEMASKS

One Size

Available with ties, earloops or headloops

Available in Sterile and Non-Sterile forms


BEARD SNOOD AND MOP CAP

One Size

Comfort fit

Product Code	Description	Size	Pack Qty	Case Qty
MNTW01100	Multi Packed Tie-On Face Mask Non-Sterile	One Size	100	10 x 100
MNTW01020	Individually Packed Tie-On Face Mask Non-Sterile	One Size	20	50 x 20
MSTW01020	Individually Packed Tie-On Face Mask Sterile	One Size	20	15 x 20
MNEW01100	Multi Packed Earloop Face Mask Non-Sterile	One Size	100	10 x 100
MNEW01020	Individually Packed Earloop Face Mask Non-Sterile	One Size	20	50 x 20
MSHW01020	Individually Packed Earloop Face Mask Sterile	One Size	50	20 x 50
MSEW01020	Individually Packed Headloop Face Mask Sterile	One Size	20	25 x 20
MNHW01100	Multi Packed Headloop Face Mask Non-Sterile	One Size	100	10 x 100
MNHW01020	Individually Packed Headloop Face Mask Non-Sterile	One Size	20	50 x 20
MNWH01ZZZ	Beard Snood White Non-Sterile	One Size	100	5 X 100
HNWH01ZZZ	Mob Cap White Non-Sterile	One Size	100	10 X 100

MATS





ACCESS MAT MEDIUM TACK NON-STERILE GUARDTACK 1

USAGE

Peel-Off access mat designed to capture contamination before entry into a controlled area

BENEFITS

Medium "Tack" formulation gives a higher level of dirt capture than GuardTack 2 whilst not inhibiting free movement

FEATURES

30 Individually numbered LDPE peel-off sheets per contamination control mat
8 mats per pack
Available in three colours (Blue / Grey / White)
Available in three sizes

Product Code	Description	Size	Pack Qty	Case Qty
ANXX08018	Contamination Control Mats Medium Tack Non-Sterile	458mm x 1145mm (18" x 45")	30 sheets/mat	8 Mats
ANXX08026	Contamination Control Mats Medium Tack Non-Sterile	660mm x 1145mm (26" x 45")	30 sheets/mat	8 Mats
ANXX08336	Contamination Control Mats Medium Tack Non-Sterile	914mm x 914mm (36" x 36")	30 sheets/mat	8 Mats

Where XX = BL for Blue, GY for Grey and WH for White



+ info see page 41



ACCESS MAT LOW TACK NON-STERILE GUARDTACK 2

USAGE

Peel-Off access mat designed to capture contamination before entry into a controlled area

BENEFITS

Low "Tack" formulation gives a lower level of dirt capture compared to GuardTack 1 but does not affect free movement

FEATURES

30 Individually numbered LDPE peel-off sheets per contamination control mat
4 mats per pack
Available in three colours (Blue / Grey / White)
Available in five sizes

Product Code	Description	Size	Pack Qty	Case Qty
ANXX04018	Contamination Control Mats Low Tack Non-Sterile	458mm x 1145mm (18" x 45")	30 sheets/mat	4 Mats
ANXX04236	Contamination Control Mats Low Tack Non-Sterile	660mm x 914mm (24" x 36")	30 sheets/mat	4 Mats
ANXX04026	Contamination Control Mats Low Tack Non-Sterile	660mm x 1145mm (26" x 45")	30 sheets/mat	4 Mats
ANXX04036	Contamination Control Mats Low Tack Non-Sterile	914mm x 1145mm (36" x 45")	30 sheets/mat	4 Mats
ANXX04336	Contamination Control Mats Low Tack Non-Sterile	914mm x 914mm (36" x 36")	30 sheets/mat	4 Mats

Where XX = BL for Blue, GY for Grey and WH for White











RECYCLING INFORMATION



TRIGGER SPRAY









Page 05 | 12

Component Name	Component Material	Disposal
 Trigger Spray Head with Dip Tube	Mixed Materials	 Non-recyclable
 500ml/950ml Bottle	HDPE	 Recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable



IPA POLYESTER POUCH WIPES









Page 09

Component Name	Component Material	Disposal
 Wipe	Polyester	 Recyclable
 Wipe Pouch	Polyester and Polyethylene blend	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable



5 LITRE RFU









Page 06 | 13

Component Name	Component Material	Disposal
 Container Cap	HDPE	 Recyclable
 5L Container	PP	 Recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable



IMS SILVER PACK POLYETHYLENE/POLYPROPYLENE WIPES











Page 10

Component Name	Component Material	Disposal
 Wipe	Polyethylene / Polypropylene	 Non-recyclable
 Wipe Pouch	OPP + PE + VMPET	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable



POLYPROPYLENE TUB WIPES









Page 07

Component Name	Component Material	Disposal
 Wipe	Polypropylene	 Recyclable
 Tub	HDPE	 Recyclable
 Tub Lid	PE	 Recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable



ISOLATOR BONNET AND TOOLS









Page 16

Component Name	Component Material	Disposal
 Bonnet	Polyester	 Recyclable
 Handle	Stainless Steel	 Recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable



POLYCELLULOSE POUCH WIPES









Page 08 | 10 | 14

Component Name	Component Material	Disposal
 Wipe	Polycellulose	 Non-recyclable
 Wipe Pouch	Polyester and Polyethylene blend	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable



DISPOSABLE MICROFIBRE MOP STERILE

Page 17 | 18 | 19

Component Name	Component Material	Disposal
 Mop	Microfibre (Polyester and Polyamide)	 Non-recyclable
 Mop Pouch (Individually Packed Only)	Polyester and Polyethylene blend	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable

**POLYCELLULOSE WIPES / ROLL WIPES**

Page 22 | 23 | 27

Component Name	Component Material	Disposal
Wipe	Hydroentangled 55% Polyester, 45% Cellulose blend	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable

**DISPOSABLE COVERALL / DISPOSABLE LAB COAT**

Page 32 | 33 | 34

Component Name	Component Material	Disposal
Coverall	White Microporous PP with PE laminate with zip and rubber elastic	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable

**POLYESTER WIPES STERILE**

Page 24 | 25 | 28

Component Name	Component Material	Disposal
Wipe	100% Polyester	Recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable

**CLEANROOM POLYESTER SOCKS**

Page 35

Component Name	Component Material	Disposal
Socks	Polyester/ HDPE	Recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable

**MICRODENIER WIPES NON-STERILE**

Page 26

Component Name	Component Material	Disposal
Wipe	100% continuous filament nylon/ polyester wipe	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable

**DISPOSABLE FACEMASKS, BEARD SNOOD AND MOB CAP**

Page 36

Component Name	Component Material	Disposal
Mobcap/Beard Snood	PP Fabric With Rubber Elastic	Non-recyclable
Facemask	Combined polymer ES (Ethylene Side By Side), Melt Blown PP (MBPP), Polyurethane, metal stripped with coating	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	

**INDIVIDUALLY PACKED GOGGLE NON-STERILE**

Page 30

Component Name	Component Material	Disposal
Goggle	Thermoplastic Rubber, Silicone, Polycarbonate	Non-recyclable
Autoclave Bag	French Paper With Blue laminated Film (PE and CPP blend)	Non-recyclable
Outer Bag	Nylon	Recyclable
Box	Cardboard	Recyclable

**ACCESS MAT**

Page 38 | 39

Component Name	Component Material	Disposal
Mat		Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable

Recycling information relates to

- UK Law and Practices
- Clean items that are free from contamination and residues



Roman Bank | Skegness | Lincolnshire | PE25 1SQ | UK
T: +44 (0)1754 767 377
E: sales@micronclean.co.uk
W: www.micronclean.com



MCB 03/2023