

# LAUNCH NEW PRODUCT

# DELTA NEUTRAL DETERGENT

Micronclean are pleased to announce the launch of a new product, Delta. Delta is a neutral detergent which has been developed for use within cleanrooms alongside our existing rotational disinfectants Alpha and Beta.

The new Annex 1 draws a clear distinction between cleaning and disinfection with a clear statement that for the disinfectant to be effective it must follow a cleaning stage.

Our Delta (Neutral Detergent) product has been developed to be fully compatible with cleanroom disinfection regimes by having the following key features:

- Low foaming formulation – this reduces the risk of streaks and staining when in use which could affect the subsequent disinfectant effectiveness.
- Neutral (pH 6.5-7.5) – to ensure that it is non-corrosive and will not damage hard surfaces.
- Supplied Pre-mixed -to ensure a homogenous distribution of the active ingredient.
- Manufactured with WFI – to ensure maximum cleanroom compatibility with cleanroom environment.
- Gamma Irradiated – to give a sterile product with a two-year shelf life.

At launch Delta will be available in the following forms: Sterile reusable laundered mops, Sterile 950ml trigger sprays and Sterile 5L RFU. In addition, we are working on other product forms which will be released in due course and will include Sterile pre-saturated wipes and Sterile disposable mops.



For more information on how Delta can help you be compliant please contact Micronclean's Cleanroom Solutions Team on: Email [SalesEnquiries@micronclean.co.uk](mailto:SalesEnquiries@micronclean.co.uk), telephone **01754 767377** or visit [www.micronclean.com](http://www.micronclean.com)

## ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
ZSND05950 - A/B/C/D/E/F/G/I	CleanGuard 5 Delta 950ml Trigger Spray Sterile	950ml	Each	15	442 x 390 x 320mm	16 Kg
ZSND05005 - A/B/C/D/E/F/G/I	CleanGuard 5 Delta 5L RFU	5L	Each	2	305 x 229 x 305mm	10 Kg