



# PRODUCT SPECIFICATION

## REUSABLE CLEANROOM SOCK

### PRODUCT INFORMATION

#### PRODUCT & SERVICE FEATURES

The Micronclean Reusable Sock service is tailored to your needs and is;

- Designed specifically for the cleanroom environment.
- Fully compliant with the requirements detailed in EU GMP Annex 1.
- Able to offer a major reduction in waste going to landfill compared with single use disposable socks.
- Fully RFID tagged so that it is compatible with our Protrack garment tracking system.

Alongside this the socks itself has not been forgotten in that it offers socks that;

- Are logo'd to support simple visual checks that show if the correct socks are being worn.
- Generate increased wearer comfort by being 'sock' shaped.
- Are available in three colour coded sizes to again generate long term wearer comfort.
- Features a synthetic Technical Fabric giving;
  - A soft feel
  - High Levels of Breathability
  - Superb Moisture Wicking

#### PRODUCT DESCRIPTION

Micronclean offer a Laundered Reuseable Cleanroom Sock service, allowing our ISO Class 4 customers amongst others to be fully compliant with EU GMP Annex 1.

#### PACKAGING

Socks are packed into single sized sealed polythene packs of up to 50 socks.

#### QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems.

#### STORAGE CONDITIONS

Store in cool, well ventilated area.

#### MATERIAL DIMENSIONS

Sock Material:

- 65% Synthetic Technical Fibre
- 22% Polyester
- 3% Rubber
- 6% Nylon
- 4% Elastic Fibre

Sock Size:

- Small - Yellow Band (up to size 6)
- Medium - Red Band (size 6 to 10)
- Large - Blue Band (size 10 to 13)