

REUSABLE CLEANROOM SOCK SERVICE

FROM MICRONCLEAN INDIA



CLEANROOM SOCK SERVICE

Micronclean offer a Laundered Reusable Cleanroom Sock Solution, allowing our ISO Class 4 customers amongst others to be compliant with EU GMP Annex 1, which states:

"Outdoor clothing (other than personal underwear) should not be brought into changing rooms leading directly to Grade B and C cleanrooms. Facility suits, covering the full length of the arms and the legs, and socks covering the feet, should be worn before entry to change rooms for Grades B and C. Facility suits and socks should not present a risk of contamination to the gowning area or processes. "




The Micronclean Reusable Sock service is tailored to your needs and is;

- Designed specifically for cleanroom environments.
- Fully compliant with the requirements detailed in EU GMP Annex 1.
- Able to offer a major reduction in waste going to landfill compared with single use disposable socks.
- Fully compatible with our Protrack garment tracking system as each sock is RFID tagged.

Alongside this the sock itself has not been forgotten in that it offers:

- Increased wearer comfort by being 'sock' shaped.
- Availability in three colour coded sizes to again generate long term wearer comfort.
- 100% synthetic fabric giving:
 - a soft feel
 - high levels of breathability
 - superb moisture wicking

CLEANROOM SOCK PROCESSING



- 01. Return of Soiled Socks** - Used socks are returned to us in their own dedicated uniquely coloured webknot bag.
- 02. Washing and Drying** - All socks are washed and tumble dried in accordance with RABC principles to ensure the thermal disinfection and cleanliness of the returned socks.
- 03. Scanning In** - All returned socks are scanned into ProTrack before being inspected.
- 04. Inspection** - The dry socks are inspected against pre-defined criteria and the replacement of non-conforming products is triggered at this point.
- 05. Packing** - Dry, inspected socks are packed into single sized sealed polythene packs of up to 50 socks, placed in the dedicated sock webknot bag ready for delivery.
- 06. Transport of Clean Socks** - The clean socks are returned as part of our delivery service alongside your existing undergarments.
- 07. Wearing** - Socks are ready for use.

CLEANROOM SOCK SERVICE

Our Cleanroom sock service is based on customer feedback and delivers sized packs of up to 50 socks alongside your existing undergarment service.

SOCK PROCESSING

Our Cleanroom sock processing is based on our existing best practice for cleanroom undergarments and delivers cleanroom socks that are processed and returned to you in a ready for immediate use.

"At Micronclean our passion is to be the first to develop new technological solutions that change the shape of the markets we serve, creating quality and efficiency for our customers."

Micronclean India has established a state-of-the-art cleanroom laundry in Bangalore. These brand new facilities have been built by Micronclean, a family-owned UK company with over 40 years' experience in the provision of contamination control.

The plant is a first in India in that it has been specifically designed to process and deliver high quality, sterile ready for immediate use cleanroom garments to support the Indian Pharma/Biotech industry.

The plant has already begun to support Pharmaceutical and Biotech producers who are seeking to export products under increased regulatory scrutiny.

Designed from the outset to incorporate world best practices to ensure compliance with the industries highest cGMP garments.



Micronclean India Pvt Ltd
Plot No. 24 & 25 | Vemagal Industrial Area | Vemagal | Kolar | 563102 | Karnataka | India
T. +91 7829111150 | E. enquiries@micronclean.in | W. www.micronclean.in

