



PRODUCT SPECIFICATION

NIPRO NEEDLE STERILE



PRODUCT INFORMATION

PRODUCT DESCRIPTION

Needles for pharmacy compounding
20G x 38mm.

PRODUCT FEATURES

Needle: stainless steel.
Cap and hub: polypropylene.
Glue: epoxy resin.

MEDICAL DEVICE CLASSIFICATION:

Not a medical device (pharmacy compounder).

PACKAGING

Primary packaging: Four layer laminate,
Easy tear notch, Alcohol and VHP
compatible.

CONNECTION TYPE:

Luer lock (also slip compatible).

TECHNICAL INFORMATION

Compatible with gamma irradiation.
Phthalate free.
Latex free.

QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016
quality management system.

STERILISATION DETAILS

Gamma Irradiation (minimum 25kGy).
Complies with applicable parts of ISO 11137
series.

SHELF LIFE

3 years from date of manufacture.

STORAGE CONDITIONS

Store in a dry, well ventilated area and
protect from direct sunlight and heat sources.

MANUFACTURING

Manufactured by Micronclean Ltd, UK

ORDERING/PACKING INFORMATION

Product Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
20NI05CC	20G Nipro Needle Sterile 25 Packs of 5	20G	5	125	285 x 160 x 200mm	0.78 kg

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please contact us.