



## PRODUCT SPECIFICATION

### DELTA 5Litre RFU STERILE

CleanGuard **5**

#### PRODUCT INFORMATION

##### PRODUCT DESCRIPTION

Sterile Neutral Detergent in RFU format.  
Available in sterile version.  
Safe to use on any surface and is colourless and odourless.  
Pre-mixed solution ensures homogenous distribution.  
Manufactured using WFI.  
pH: 6.5 - 7.5

##### PRODUCT FEATURES

Hard surface cleaning agent.  
5 Litre container complete with a screw on air tight cap.  
Convenient carrying size and weight.  
Visible liquid level.  
Easy to use.

#### TECHNICAL INFORMATION

##### MANUFACTURING

The product is filled and capped in a Grade B (in use) cleanroom.  
Packaging is performed in a Grade C (in use) cleanroom.  
The disinfectant is filtered through a 0.2 micron filter at the point of fill.

##### PACKAGING

Primary packaging: 5 Litre container complete with a screw on air tight cap.  
Secondary packaging: double bagged in polythene bags.  
2 double-bagged bottles in the box.

##### STERILISATION DETAILS

Gamma Irradiation (minimum 25kGy).

##### QUALITY ASSURANCE

Product supplied under BS EN ISO 9001:2015 registered systems.

##### SHELF LIFE

2 years from date of manufacture.

##### STORAGE CONDITIONS

Store in a dry, well-ventilated area in ambient conditions.

#### ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
	ZSND05005	Delta 5L RFU Sterile	5L	Each	2 x 5L	290 x 215 x 290mm	10 Kg

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please contact us.