



## PRODUCT SPECIFICATION

### STAINLESS STEEL ISOLATOR CLEANING TOOL NON-STERILE

PureGuard **7**

#### PRODUCT INFORMATION

##### PRODUCT DESCRIPTION

Lightweight isolator tool with 360° rotation designed specifically to help clean inside of isolators and LAF cabinets.

##### PRODUCT FEATURES

Designed for use in conjunction with PureGuard 6 Dry Dual Layer Polyester Isolator Bonnet Sterile.

Two handle lengths are available 35cm (Short) and 61cm (Long).

##### PACKAGING

Sold Individually.

##### STERILISATION DETAILS

Supplied Non-sterile.

Suitable for gamma irradiation, VHP sterilisation and autoclaving.

##### QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems.

##### DIMENSIONS

Overall width = 10cm

Handle length = 35cm or 61cm

##### MATERIALS

Base plate material: T304 stainless steel

Swivel wire material: T304 stainless steel

Tube material:

6063A T6 anodised AA10 Aluminium

Stainless steel tube plug material:

T304 stainless steel

Plastic tube bung material: Polyamide

#### ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Case Qty
ZNSS01035		Stainless Steel Short Handle Isolator Cleaning Tool Non-Sterile	35 cm	EA	EA
ZNSS01061		Stainless Steel Long Handle Isolator Cleaning Tool Non-Sterile	61 cm	EA	EA
	ISED014	Dry Dual Layer Polyester Isolator Bonnet Sterile	One Size	Each	40