



# PRODUCT SPECIFICATION

## MINI-SPIKE 2 DISPENSING PINS

### PRODUCT INFORMATION

#### PRODUCT DESCRIPTION

A vented dispensing pin for injecting / withdrawing fluids from single and multidose containers. The large-sized 0.2 µm bacteria retentive air filter protects against microbial contamination and reduces the risk of exposure with toxic aerosols.

#### PRODUCT FEATURES

Grip plate with a female luer connector for a stable connection between container and syringe.

Female luer connector is featured with a two-way swabable valve to ensure clean, leak-proof and closed drug admixture.

Micro spike with shorter length to minimize residual volumes in small vials.

#### Materials

Grip plate: Styrene-Acrylonitrile.

Air filter: Acrylic Copolymer on nonwoven nylon support.

Piercing spike: Styrene - Acrylonitrile / Acrylonitrilebutadiene-styrene plus silicone lubricant.

Spike protection: Polyethylene (plus masterbatch green).

Swabable valve: Polycarbonate (plus masterbatch red).

#### PACKAGING

Primary packaging: Four layer laminate, complete with easy tear notch, Alcohol and VHP compatible.

#### TECHNICAL INFORMATION

Compatible with gamma irradiation.

BSE/TSE - considered not to present any risk.\*

Phthalate free.

Latex free.

PVC free.

\*This product does not contain tissues or cells of animal origin or their derivatives as referred to in Regulation (EU) 722/2012. Tallow derived products used in this product fulfill the requirements laid down in Regulations 1069/2009/EC and 142/2011/EC and the 'Note for Guidance EMEA/410/01, rev 3.

#### QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016 quality management system.

#### STERILISATION DETAILS

Gamma irradiation min. 25kGy.

Complies with applicable parts of ISO 11137 series.

#### SHELF LIFE

3 years from date of manufacture.

#### STORAGE CONDITIONS

Store in a dry, well ventilated area and protect from direct sunlight and heat sources.

#### MANUFACTURING

Micronclean Ltd, Holly Road, Skegness, Lincolnshire, UK

## ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
CDPM25CC	Mini-Spike 2 Dispensing Pin Sterile Pack of 25	n/a	25	100	330 x 280 x 400mm	3.10 KG

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please contact us.