



PRODUCT SPECIFICATION

CHANGEABLE GOGGLE LENS NON-STERILE

VisionGuard **1**

PRODUCT INFORMATION

PRODUCT DESCRIPTION

The ONACL0110 Lens is the changeable goggle lens for the ONACL0112 Autoclavable goggle.
Supplied Non-Sterile.

PRODUCT FEATURES

Polycarbonate anti-fog lens.
Fully compatible with ONACL0112 Autoclavable goggle.
Supplied Non - Sterile individually packed in an outer polyethylene pouch.
Autoclave life up to 40 cycles.
Can be autoclaved on its side or when installed in a goggle housing.

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001:2015 registered systems.

PHYSICAL DIMENSIONS

Lens Width: 80mm
Lens Length: 178mm
Lens Thickness : 2mm
Weight (per lens): 0.035 kg

STERILISATION (IF APPLICABLE)

Supplied in a Non-Sterile condition but can be autoclaved.

SHELF LIFE

2 years Non-Sterilised.

STORAGE CONDITIONS

Store in ambient conditions.

RECYCLING INFORMATION

Component Name:	Component Material:	Disposal:	Guidance:
Lens	Polycarbonate	Non-recyclable	Dispose of in line with local guidelines
Pouch	LDPE	Recyclable	
Cardboard Box	Cardboard	Recyclable	

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
ONACL0110	Changeable Goggle Lens Non-Sterile	One Size	Each	10 packs of 1	213 x 204 x 100mm	0.50 kg