


# REUSABLE CLEANROOM SOCK SERVICE

---



## CLEANROOM SOCK SERVICE



Micronclean offer a Laundered Reusable Cleanroom Sock Solution, allowing all our customers to be compliant with EU GMP Annex 1, which states:

*"Outdoor clothing (other than personal underwear) should not be brought into changing rooms leading directly to Grade B and C cleanrooms. Facility suits, covering the full length of the arms and the legs, and socks covering the feet, should be worn before entry to change rooms for Grades B and C. Facility suits and socks should not present a risk of contamination to the gowning area or processes. "*



**The Micronclean Reusable Sock service is tailored to your needs and is;**

- Designed specifically for cleanroom environments.
- Fully compliant with the requirements detailed in EU GMP Annex 1.
- Able to offer a major reduction in waste going to landfill compared with single use disposable socks.
- Fully compatible with our Protrack garment tracking system as each sock is RFID tagged.


**Alongside this the sock itself has not been forgotten in that it offers:**

- Micronclean logos to support simple visual checks that show if the correct socks are being worn.
- Increased wearer comfort by being:
  - 'Sock' shaped
  - Available in three colour coded sizes
- A polyester blend fabric giving:
  - A soft feel
  - High levels of breathability
  - Superb moisture wicking

---

## CLEANROOM SOCK PROCESSING

---



- 01. Return of Soiled Socks** - Used socks are returned to us in their own dedicated uniquely coloured webknot bag.
- 02. Washing and Drying** - A validated cleanroom laundering process where the socks are washed and tumble dried in a dedicated cleanroom facility.
- 03. Scanning In** - All returned socks are scanned into ProTrack before being inspected.
- 04. Inspection** - The dry socks are inspected against pre-defined criteria and the replacement of non-conforming products is triggered at this point.
- 05. Packing** - Dry, inspected socks are packed into single sized sealed polythene packs of up to 50 socks, placed in the dedicated sock webknot bag ready for delivery.
- 06. Transport of Clean Socks** - The clean socks are returned as part of our delivery service alongside your existing undergarments.
- 07. Wearing** - Socks are ready for use.

### BECAUSE EVERY STEP MATTERS

Driven by customer feedback, our Cleanroom Sock Service enhances our industry-leading garment laundry program, delivering a complete, head-to-toe solution for controlled environments. Each sock is processed in a dedicated cleanroom laundry facility using a validated wash/dry cycle that is proven to achieve a 5-log reduction in bioburden. This ensures every sock meets the highest standards of contamination control. In fact, these might be the cleanest socks you'll ever wear. [[Read the Reusable Cleanroom Sock Report](#)]




**“Our aim is to be the most sustainable specialist textile and cleanroom consumables business on the planet.”**

Micronclean is committed to contributing to a sustainable future based on the three pillars of:

- Protecting our climate
- Creating our circular economy
- Delivering a positive impact for our people and communities

**microngreen**  
our environmental future



ISSUE 3 | DATE 03/06/2025



Roman Bank | Skegness | Lincolnshire | PE25 1SQ | UK  
T: +44 (0)1754 767 377  
E: sales@micronclean.co.uk  
W: www.micronclean.com

