




STEWARDSHIP / KNOWLEDGE / INNOVATION / **EXCELLENCE**

The Micronclean SKIEs initiative - at the very heart of everything we do at Micronclean



MICRONCLEAN STERILE SYRINGE PACKS

A comprehensive range of sterile multi-packed syringes.
100% quality checked and packed in our new
MicronDevices cleanroom



MICRONCLEAN STERILE SYRINGE PACKS

PRODUCT DESCRIPTION

MICRONCLEAN PACK THE SYRINGES
IN OUR OWN MICRONDEVICES
(ISO CLASS 6) CLEANROOM IN
SKEGNESS, UK



KEY FEATURES AND BENEFITS

| FEATURES | BENEFITS |
|--|---|
| Single batches of syringes | Simplified paperwork with increased efficiency |
| Consolidated Certificate of Conformity | Single document, easy to file and cross reference |
| Quad laminate pouch | Transparent, strong, easy to open, works for VHP and alcohol transfer process |
| 100% tipped syringes | Secure in usage |
| Broad heat seal | Burst resistant device packaging |
| Large stocks held by Micronclean | Security of supply |
| Logistics by Micronclean | Easy to set service up with low or no delivery charges |
| Storage | Validated shelf life |
| Packed by Micronclean | Simplified supply chain with increased efficiency |

QUALIFICATION

Micronclean are registered to BS EN ISO 13485:2003/12 which covers every aspect of sourcing, manufacturing, sterilisation and distribution of the sterile pack range. Micronclean hold a detailed technical file covering all products in the range.

TRACEABILITY AND STERILISATION

Full batch traceability and product certification is as standard. As each pouch of products has a single traceable batch number (traceable to all the sub-components batch numbers) the work of the technician in transcribing information is significantly reduced as is the potential for error. Sterilisation is by a fully validated gamma irradiation process. The product is supplied with a consolidated certificate of conformity.

PACKAGING

The syringe and component packs are supplied in transparent, easy-tear, quad laminate pouches which are VHP and alcohol impermeable, making them suitable for most transfer processes. In addition the syringes are all supplied with a 'tip protector' to protect the critical tip during storage and transit.

The syringe packs are double bagged to ensure an easy and rapid transfer cycle.

ORDERING & DELIVERY

Easy order placement and stock holding is managed by Micronclean's powerful stock inventory software. Stocking and delivery is direct from Micronclean's state-of-the-art storage facilities and own fleet of delivery vehicles enabling us to deliver right to the point of requirement.

CUSTOM PACKS

Micronclean are happy to discuss any requirements you may have for custom made packs which incorporate other sterile ancillaries.



MICRONCLEAN < STERILE SYRINGE PACKS

ORDERING INFORMATION

| BD Syringe Size Use Volume | Product Code | Description | Items per pouch | Pouches per carton |
|-------------------------------|--------------|--|--------------------|-----------------------|
| 1ml | 01ML05CC | 1ml luer lock syringe plus tip protector | 5 | 20 |
| | 01ML10CC | 1ml luer lock syringe plus tip protector | 10 | 10 |
| | 01ML25CC | 1ml luer lock syringe plus tip protector | 25 | 10 |
| | 01ML25SC | 1ml luer lock syringe plus barcode and caps sleeved | 25 | 12 |
| 3ml | 03ML05CC | 3ml luer lock syringe plus tip protector | 5 | 20 |
| | 03ML10CC | 3ml luer lock syringe plus tip protector | 10 | 10 |
| | 03ML25CC | 3ml luer lock syringe plus tip protector | 25 | 10 |
| | 03ML25CH | 3ml luer lock syringe plus tip protector and cap | 25 | 10 |
| | 03ML50CC | 3ml luer lock syringe plus tip protector | 50 | 10 |
| 5ml | 05ML25SC | 3ml luer lock syringe plus barcode and caps sleeved | 25 | 12 |
| | 05ML05CC | 5ml luer lock syringe plus tip protector | 5 | 20 |
| | 05ML10CC | 5ml luer lock syringe plus tip protector | 10 | 10 |
| | 05ML25CC | 5ml luer lock syringe plus tip protector | 25 | 10 |
| | 05ML25CH | 5ml luer lock syringe plus tip protector and cap | 25 | 10 |
| 10ml | 05ML50CC | 5ml luer lock syringe plus tip protector | 50 | 10 |
| | 10ML05CC | 10ml luer lock syringe plus tip protector | 5 | 10 |
| | 10ML10CC | 10ml luer lock syringe plus tip protector | 10 | 10 |
| | 10ML25CC | 10ml luer lock syringe plus tip protector | 25 | 10 |
| | 10ML25CH | 10ml luer lock syringe plus tip protector | 25 | 10 |
| 20ml | 10ML25SC | 10ml luer lock syringe plus tip protector and cap | 25 | 8 |
| | 10ML50CC | 10ml luer lock syringe plus barcode and caps sleeved | 50 | 5 |
| | 20ML05CC | 20ml luer lock syringe plus tip protector | 5 | 10 |
| | 20ML10CC | 20ml luer lock syringe plus tip protector | 10 | 10 |
| | 20ML25CC | 20ml luer lock syringe plus tip protector | 25 | 10 |
| 30ml | 20ML25SC | 20ml luer lock syringe plus barcode and caps sleeved | 25 | 8 |
| | 20ML50CC | 20ml luer lock syringe plus tip protector | 50 | 10 |
| | 30ML05CC | 30ml luer lock syringe plus tip protector | 5 | 10 |
| | 30ML10CC | 30ml luer lock syringe plus tip protector | 10 | 10 |
| | 30ML25CC | 30ml luer lock syringe plus tip protector | 25 | 10 |
| 50ml | 30ML25CH | 30ml luer lock syringe plus tip protector and caps | 25 | 10 |
| | 30ML25SC | 30ml luer lock syringe plus barcode and caps sleeved | 25 | 8 |
| | 50ML05CC | 50ml luer lock syringe plus tip protector | 5 | 10 |
| | 50ML10CC | 50ml luer lock syringe plus tip protector | 10 | 10 |
| | 50ML25CC | 50ml luer lock syringe plus tip protector | 25 | 8 |
| | 50ML25CH | 50ml luer lock syringe plus tip protector and caps | 25 | 8 |
| | 50ML25SC | 50ml luer lock syringe plus barcode and caps sleeved | 10 | 6 |
| | 50ML300C | 50ml luer lock syringe plus barcode and caps sleeved | 25 | 12 |

Ancillary Items - Qualified for VHP and Alcohol Transfer Processes (1)

| Double Bagged Product | Order Code | No. per pouch | Pouches per box |
|-----------------------------|------------|------------------|--------------------|
| Baxa luer lock syringe caps | CAPS10CC | 10 | 10 |
| Baxa luer lock syringe caps | CAPS50CC | 50 | 10 |
| Chemo dispensing pin | CDPB25CC | 25 | 4 |
| Large volume dispensing pin | LVDP10CC | 10 | 10 |
| B Braun needle 19g x 1.5" | 19GN50CC | 50 | 10 |
| B Braun needle 21g x 1.5" | 21GN50CC | 50 | 10 |
| Terumo needle 19g x 1.5" | 19TR50CC | 50 | 10 |
| Terumo needle 21g x 1.5" | 21TR25CC | 25 | 20 |

| Double Bagged Product | Order Code | No. per pouch | Pouches per box |
|--------------------------|------------|------------------|--------------------|
| Terumo needle 21g x 1.5" | 21TR50CC | 50 | 10 |
| Filter needle | 19FN25CC | 25 | 10 |
| Codan dispensing pins | COSP25CC | 25 | 10 |
| Chemo pack | CHTH01CC | N/A | 8 |
| Barcode labels | BCL200CC | 200 | 10 |
| Environmental labels | ENV200CC | 200 | 20 |
| Pens | PEN005CC | 5 | 20 |
| Waste bags | SBAG10CC | 10 | 10 |

(1) Other sanitisation methods may be appropriate if qualified by the end user



STEWARDSHIP / KNOWLEDGE / **INNOVATION** / EXCELLENCE

Supporting cGMP Compliance Across the Globe / Exporting Excellence Worldwide



micronclean
RD&T
RESEARCH DEVELOPMENT & TESTING UNIT

PRODUCT DEVELOPMENT

This brochure contains products stocked for immediate despatch from the UK warehouse. However, Micronclean also operate an RD&T Unit that is focussed on creating innovative new products and solutions to serve our customers worldwide. If there is a product that you require but cannot see it in this brochure, please call us today, we will be more than happy to put the team to work finding or creating a bespoke product or solution that meets the needs of your application and environment

micronclean
QUALITY TODAY - INNOVATION TOMORROW

Roman Bank / Skegness / Lincolnshire / PE25 1SQ

T: +44 (0)1754 767 377

E: enquiries@micronclean.co.uk

W: www.micronclean.com



For more information about our SKIEs initiative, please ask for a copy of our SKIEs brochure.