



FAQs

UPGRADED WSCW97ZZZ VERIGUARD 9 IMS SILVER PACK WIPE

Q1. Why are you changing the specification of the ‘Silver Pack’ pre-saturated wipe?

As part of our ongoing commitment to you as a customer we are always looking to improve our products to improve your efficiency. As a result, we are upgrading the wipe material and improving the way the product is folded in the silver pack pouch.

Q2. What is the timescale for the introduction of the improved “Silver Pack” pre-saturated wipe?

We had planned to upgrade this product in early 2020 but due to unforeseen circumstances these improved wipes have been delayed until June 2020.

To cover a potential shortfall in stock of the existing Silver Pack wipe we have re-introduced WSVG07020 VeriGuard 1 IMS Polypropylene Tub Wipe Sterile. The product specification and trial samples of this product are available.

Q3. Why have you upgraded the wipe material?

By moving from a Polypropylene/Viscose fibre to a Polyethylene/Polypropylene fibre for the same weight of wipe we have improved its absorbency, giving you the ability to wipe a larger area as the wipe does not dry out as quickly.

Note any changes to your cleanroom regime will have to be locally verified against your existing operating procedures.

Q4. Why have you changed the wipe fold?

By replacing the interleaved Z-fold with a non-interleaved Z-fold we have made it easy for an operator to remove a single wipe from the pouch as well as it making it easier to reseal the pack as the opening will not be obstructed by the next wipe.

Q5. Have you made any other changes to the wipe?

No other changes to the wipe or its processing so that the new product retains:

- The existing easy peel opening
- The resealable pouch for improved alcohol retention
- The same grade of IMS and
- Irradiation in the UK

Q6. With all the changes why have you retained the same Product code and Family Name and Number as the original product?

We have retained the same product code and family name and number as the original product as we see this as helping all our customers make the transition from the old to the new improved product. However, we have added a regional code suffix to the product code to help with BPR registration requirements going forward. In addition we have taken the opportunity to rename the product in line with our new naming conventions (see below for details).

Product Code	Product Family and Number	Product Name
WSCW97ZZZ -A/1	VeriGuard 9	IMS Silver Pack Polyethylene/Polypropylene Z-folded Wipes Sterile