

MICRONCLEAN COMPLIANCE 100™ WIPES SPECIFICATION

Product Code :	WSPW60300
Product Description :	Sterile cleanroom wipes made from twin layer continuous filament polyester with sealed edges
Pack Size :	5 wipes / pack
Pouch Format :	Micro25

Technical Specification

Fabric :	100% polyester; twin-layer construction; 264 g/m ² weight (tbc); tested to ASTM-D-629
Wipe Size :	300mm x 300mm
Decontamination :	Aqueous washing incorporating thermal disinfection to minimum 71°C / 3 minutes; product packed in ISO 14644-1 Class 4 (at rest), EU GMP Grade B cleanroom
Sterilisation :	Validated sterilisation by gamma irradiation (minimum 25 kGy)
Shelf-Life :	Two years from manufacturing date
Quality Assurance :	Product manufactured under BS EN ISO 9001:2000 & BS EN ISO 13485:2003 registered systems