

# CONSUMABLES

UK BROCHURE :: 2025





## CONTENTS

INTRODUCTION	01	DISPOSABLE LABCOAT NON-STERILE	34
REGIONAL ORDERING AND LABELS	02	SINGLE USE CLEANROOM COMPATIBLE SOCK NON-STERILE	35
<b>ALCOHOLS</b>	<b>03</b>	DISPOSABLE FACEMASKS, BEARD SNOOD AND MOB CAP NON-STERILE	36
IPA HAND SANITISER	04	<b>GLOVES</b>	<b>37</b>
IPA AND IMS TRIGGER SPRAY	05	GLOVE DECISION TREE	38
IPA AND IMS RFU 5L	06	POLYCHLOROPRENE 300mm HAND SPECIFIC GLOVE STERILE	39
IPA POLYPROPYLENE TUB WIPE	07	LATEX 300mm HAND SPECIFIC GLOVE STERILE	40
IPA AND IMS POLYCELLULOSE POUCH WIPE	08	NITRILE 300mm HAND SPECIFIC GLOVE STERILE	41
IPA POLYESTER POUCH WIPE	09	NITRILE 405mm HAND SPECIFIC GLOVE STERILE	42
IMS SILVER PACK POLYETHYLENE/POLYPROPYLENE PRE-SATURATED WIPE	10	NITRILE 230MM AMBIDEXTROUS GLOVE BLUE BOXED NON-STERILE	43
<b>DISINFECTANTS AND DETERGENT</b>	<b>11</b>	LATEX 300MM AMBIDEXTROUS GLOVE BULK PACKED NON-STERILE	44
ALPHA, BETA AND DELTA TRIGGER SPRAY	12	NITRILE 300MM AMBIDEXTROUS GLOVE BULK PACKED NON-STERILE	45
ALPHA, BETA AND DELTA AND DI 5L RFU	13	NITRILE 300mm AMBIDEXTROUS GLOVE BULK PACKED NON-STERILE	46
ALPHA, BETA AND DELTA POLYCELLULOSE POUCH WIPE	14	<b>STATIONERY</b>	<b>47</b>
<b>DISPOSABLE MOP SYSTEM</b>	<b>15</b>	CLEANROOM BALLPOINT PEN AND MARKER PEN BLACK STERILE	48
ISOLATOR BONNET AND TOOL	16	CLEANROOM PAPER NON-STERILE	49
IPA AND IMS DISPOSABLE MICROFIBRE MOP STERILE	17	CLEANROOM MUNISING PAPER STERILE AND NON-STERILE	50
ALPHA, BETA AND DELTA DISPOSABLE MICROFIBRE MOP STERILE	18	<b>OVERSHOES</b>	<b>51</b>
DRY DISPOSABLE MICROFIBRE MOP STERILE	19	POLYTHENE OVERSHOE BLUE NON-STERILE	52
WASTE BAGS, MOP HEADS, TROLLEY AND OTHER ACCESSORIES	20	PVC OVERSHOE NON-STERILE	53
<b>DRY WIPES</b>	<b>21</b>	OVERSHOE WITH TEXTURED PVC SOLE AND ANTI-SLIP NON-STERILE	54
WIPE DECISION TREE	22	<b>MATS</b>	<b>55</b>
POLYCELLULOSE WIPE STERILE	23	ACCESS MAT MEDIUM TACK NON-STERILE	56
POLYCELLULOSE WIPE NON-STERILE	24	ACCESS MAT LOW TACK NON-STERILE	57
POLYESTER WIPE STERILE	25	BLUE WALK ACCESS MAT WASHABLE NON-STERILE	58
POLYESTER WIPE NON-STERILE	26	GREY PRO ACCESS MAT WASHABLE NON-STERILE	59
MICRODENIER WIPE NON-STERILE	27	BLUE SELF-INSTALL ACCESS MAT NON-STERILE	60
POLYCELLULOSE ROLL WIPE NON-STERILE	28	<b>MICRONGREEN</b>	<b>61</b>
<b>GOGGLES</b>	<b>29</b>	<b>RECYCLING INFORMATION</b>	<b>63</b>
INDIVIDUALLY PACKED GOGGLE NON-STERILE	30	<b>MICRONCLEAN ASSEMBLY</b>	<b>67</b>
<b>CLOTHING</b>	<b>31</b>		
DISPOSABLE COVERALL WITH INTEGRAL HOOD	32		
DISPOSABLE COVERALL WITH INTEGRAL HOOD AND FEET STERILE	33		



---

## INTRODUCTION

---

*"At Micronclean, it is our passion to be the first to develop new technological solutions that transform the nature of the markets we serve, creating quality and efficiency for our customers."*

For over 40 years, Micronclean has focused on contamination control products for companies in the pharmaceutical, healthcare, microelectronics, aerospace, and semiconductor industries.

Therefore, we are incredibly proud of our range of cleanroom consumables, which we have sold in over 50 countries worldwide over the past 12 years.

Much of our range is manufactured in the UK under cleanroom conditions. For the products we don't manufacture, we have full control over the supply chain and set the specifications. Our manufacturing capability means we can offer customized products and maintain sufficient stock levels to meet customer demand. Our customers can also audit our facilities.

### **MicronGreen: Our Commitment to Sustainability**

At Micronclean, sustainability is a core value that guides everything we do. Through our MicronGreen program, we pursue the vision of *"becoming the most sustainable specialist textiles and cleanroom consumables company in the world"*.

Since its launch in 2022, MicronGreen has made significant progress in reducing our environmental footprint. Between 2021 and 2023, we reduced our Scope 1 and Scope 2 emissions by 20% by implementing energy efficiency measures, installing solar panels, and fully transitioning to 100% renewable energy.

Our commitment to waste reduction has led to innovative packaging solutions. By combining two separate labels into a single wipe-top closure, we were able to reduce our annual plastic consumption by 200 kg. Furthermore, switching from white-branded cartons to eco-friendly brown cartons with recycled fiber content resulted in a 32% reduction in the use of virgin cardboard.

**Angela Holland**  
Group Product Manager Consumables

For more information on any of our product lines please visit our website [www.micronclean.com](http://www.micronclean.com) or contact us directly on [cs@micronclean.co.uk](mailto:cs@micronclean.co.uk) or +44 (0) 1754 767377

# REGIONAL ORDERING AND LABELS

## MICRONCLEAN UPDATES ITS PRODUCT LABELS TO DEAL WITH DEMAND FROM EXPORT MARKETS.

Following increasing demand for our market leading disinfectant and cleaning products, we have had to enhance our existing product labels to ensure that we remain compliant.

As a result, we have generated an enhanced label which enables us to remain compliant with labeling regulations and allows us to carry up to three languages in addition to English.

So, the result of these changes is that the product code for our existing CleanGuard 1 IPA 500ml Trigger Spray Sterile will migrate from the current ZSAL64499 to ZSAL64499-A for region A orders, ZSAL64499-B for region B orders and so on. Note this regionalisation only affects products with active chemicals.

The new look labels with the associated regional codes for products will be rolled out progressively across our whole consumable range.

### Regional Labels in Use

For a UK consumable customer, the default option will always be that they order Region A products and they will normally receive Region A products. However, as the English language is on the majority of labels and labelling regulations are common across regions A-E, J and H, and we might on occasion substitute the same product with a region K label for example.

For a Micronclean customer in Italy for instance, the product should always be ordered for Region B and to comply with the local language regulations we will supply the product with Region B labelling.

For regions F, G, H and I, the labelling regulations are country specific and as such no substitutions are possible.

Please note if you have any questions please do not hesitate to contact us.

REGIONS	A	UK   Ireland   Germany   Austria
	B	UK   Spain   Portugal   Italy
	C	UK   France   Belgium   Holland
	D	UK   Sweden   Denmark   Norway
	E	UK   Greece   Poland
	F	China
	G	Australia
	H	Brazil
	I	India
	J	UK   Croatia   Turkey
	K	UK   Slovenia   Hungary

# ALCOHOLS





+ info see page 64



## HAND SANITISER

### CLEANGUARD 1

#### USAGE

Designed as a sanitising hand rub for use in a professional environment

#### BENEFITS

70% alcohol content  
 Manufactured with WFI  
 UK BPR Registered  
 Validated Performance:  
 - EN 13727:2012+A2:2015 – Pass  
 - EN 13624:2013 – Pass

#### FEATURES

Transparent bottle  
 Simple push down dispensing pump  
 Fragrance-free  
 Non-sticky  
 Quick drying  
 Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
ZNIP50025	IPA 500ml Hand Sanitiser Non-Sterile	500ml	Each	24 x 500ml

## RECYCLING INFORMATION



+ info see page 64



## ALCOHOL TRIGGER SPRAY

### CLEANGUARD 1 AND 2

#### USAGE

Designed for use in Cleanrooms  
Cleanroom transfer process  
Disinfection of hard surfaces  
Use within Isolators / LAFs

#### BENEFITS


No bag in the bottle  
70% alcohol blended with WFI quality water  
Filtered to 0.2µm at point of fill  
Mechanically sealed trigger head  
Low centre of gravity bottle design  
Transparent bottle  
Long trigger handle, ease of use  
12 week sterility from first use

#### FEATURES

Trigger spray with spray / off feature  
2 year shelf life  
Gamma irradiation stable  
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
ZSAL64499	IPA 500ml Trigger Spray Sterile	500ml	Each	24 x 500ml
ZSAL64950	IPA 950ml Trigger Spray Sterile	950ml	Each	15 x 950ml
ZNAL64950	IPA 950ml Trigger Spray Non-Sterile	950ml	Each	15 x 950ml
ZSAL65499	IMS 500ml Trigger Spray Sterile	500ml	Each	24 x 500ml
ZSAL65950	IMS 950ml Trigger Spray Sterile	950ml	Each	15 x 950ml
ZNAL65950	IMS 950ml Trigger Spray Non-Sterile	950ml	Each	15 x 950ml

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



RECYCLING INFORMATION

+ info see page 64



# IPA AND IMS 5L RFU CLEANGUARD 1 AND 2

**USAGE**

Disinfection of hard surfaces

**BENEFITS**

- 70% Alcohol blended with WFI quality water
- Bulk container for larger surface areas
- Easy to use with bucket mopping systems
- Filtered to 0.2µm at point of fill

**FEATURES**

- Tamper evident PP vented cap
- Gamma irradiation stable
- 2 year shelf life
- Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
ZSAL64005	IPA 5 Litre RFU Sterile	5L	Each	2 x 5L
ZNAL64005	IPA 5 Litre RFU Non-Sterile	5L	Each	2 x 5L
ZSAL65005	IMS 5 Litre RFU Sterile	5L	Each	2 x 5L

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



## RECYCLING INFORMATION



+ info see page 64



## IPA POLYPROPYLENE TUB WIPE

### VERIGUARD 1

#### USAGE

Hard surface wiping particularly suited to outer areas

#### BENEFITS

Designed for low grade cleanroom areas  
Designed for high volume usage  
Long lasting alcohol retention  
Re-fillable pouch available

#### FEATURES

Continuous filament thermobonded nonwoven polypropylene wipe  
18gsm wipe weight  
2 years shelf life  
Sterile version only supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
WSVG05020	IPA Polypropylene Tub Wipe Sterile	200mm x 200mm	200	12 x 200
WNVG05020	IPA Polypropylene Tub Wipe Non-Sterile	200mm x 200mm	200	12 x 200

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

RECYCLING INFORMATION

+ info see page 64



# IPA AND IMS POLYCELLULOSE POUCH WIPE

## VERIGUARD 2

### USAGE

Hard surface wiping particularly suited to transfer processes

### BENEFITS

- Increased alcohol retention from re-sealable pouch
- High absorbency
- Strong multifunctional wipe
- Small pack size available for small batch manufacturing, large pack size suitable for large batch manufacturing

### FEATURES

- White hydro-entangled polycellulose wipe
- 68gsm wipe weight
- 2 years shelf life
- Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
WSVG04023	IPA Polycellulose C-folded Pouch Wipes Sterile	230mm x 230mm	100	8 Packs of 100
WSCW92009	IPA Sessional Polycellulose C-folded Pouch Wipes Sterile	230mm x 230mm	15	40 Packs of 15
WNVG04023	IPA Polycellulose C-folded Pouch Wipes Non-Sterile	230mm x 230mm	100	8 Packs of 100
WSVG09023	IMS Polycellulose C-folded Pouch Wipes Sterile	230mm x 230mm	100	8 Packs of 100
WSCW91009	IMS Sessional Polycellulose C-folded Pouch Wipes Sterile	230mm x 230mm	15	40 Packs of 15

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

RECYCLING INFORMATION



+ info see page 64



## IPA POLYESTER POUCH WIPE VERIGUARD 3

### USAGE

Cleanroom hard surface wipe  
Use in Isolators / Laminar flow cabinets

### BENEFITS

Increased alcohol retention from re-sealable pouch  
Laser cut edges  
Strong multifunctional wipe  
Small pack size available for small batch manufacturing, medium pack size suitable for larger batch manufacturing

### FEATURES

100% continuous filament polyester wipe  
2 years shelf life  
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
WSVG06030	IPA Sessional Polyester V-folded Pouch Wipe Sterile	300mm x 300mm	15	30 Packs of 15
MCZSLA50S	IPA Polyester Laundered Non-folded Pouch Wipes Sterile	300mm x 300mm	50	5 Packs of 50
WNVG06030	IPA Sessional Polyester V-folded Pouch Wipe Non-Sterile	300mm x 300mm	15	30 Packs of 15
MCZNLA50N	IPA Polyester Laundered Non-folded Pouch Wipes Non-Sterile	300mm x 300mm	50	5 Packs of 50

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



RECYCLING INFORMATION

+ info see page 64



# IMS SILVER PACK POLYETHYLENE / POLYPROPYLENE PRE-SATURATED WIPE VERIGUARD 9

**USAGE**

Cleanroom hard surface wipe

**BENEFITS**

- Increased alcohol retention from re-sealable pouch
- Designed for large batch manufacturing
- Blended with de-ionised water (DI)
- Long lasting alcohol retention

**FEATURES**

- Nonwoven interleaved wipes
- 2 years shelf life
- Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
WSCW97ZZZ	IMS Silver Pack Polyethylene / Polypropylene Z-folded Wipes Sterile	200mm x 220mm	200	8 Packs of 200

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

# DISINFECTANTS





+ info see page 64



## ALPHA, BETA AND DELTA TRIGGER SPRAY

### CLEANGUARD 3, 4 AND 5

#### USAGE

Disinfection and cleaning of hard surfaces

Annex 1 states that for disinfection to be effective it must be a two-stage process. The first stage should be a process to remove surface contamination and the second stage should be disinfection using a rotational system based on two types of disinfectants which have different modes of action

Micronclean has developed Delta as a sterile neutral detergent to remove surface contamination. This should be followed by Alpha and Beta as sterile rotational disinfectants which have different modes of action and prevent the long-term build-up of resistance

#### BENEFITS

- No bag in the bottle
- Blended with WFI quality water
- Filtered to 0.2µm at point of fill
- Mechanically sealed trigger head
- Low centre of gravity bottle design
- Transparent bottle
- Long trigger handle, ease of use

#### FEATURES

- Trigger spray with spray / off feature
- 2 year shelf life
- Gamma irradiation stable
- Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
ZSAP01500	Alpha 500ml Trigger Spray Sterile	500ml	Each	6 x 500ml
ZSBT02500	Beta 500ml Trigger Spray Sterile	500ml	Each	6 x 500ml
ZSND05950	Delta 950ml Trigger Spray Sterile	950ml	Each	15 x 950ml

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

## RECYCLING INFORMATION



+ info see page 64



## ALPHA, BETA, DELTA AND DI 5L RFU CLEANGUARD 3, 4, 5 AND 6

### USAGE

Disinfection and cleaning of hard surfaces and floors

Annex 1 states that for disinfection to be effective it must be a two-stage process. The first stage should be a process to remove surface contamination and the second stage should be disinfection using a rotational system based on two types of disinfectants which have different modes of action

Micronclean has developed Delta as a sterile neutral detergent to remove surface contamination. This should be followed by Alpha and Beta as sterile rotational disinfectants which have different modes of action and prevent the long-term build-up of resistance

### BENEFITS

Blended with WFI quality water  
Bulk container for larger surface areas  
Easy to use with bucket mopping systems  
Filtered to 0.2µm at point of fill

### FEATURES

Tamper evident PP vented cap  
Gamma irradiation stable  
2 year shelf life  
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
ZSAP03005	Alpha 5 Litre RFU Sterile	5L	Each	2 x 5L
ZSBT04005	Beta 5 Litre RFU Sterile	5L	Each	2 x 5L
ZSND05005	Delta 5 Litre RFU Sterile	5L	Each	2 x 5L
ZSDI82005	DI Water 5 Litre RFU Sterile	5L	Each	2 x 5L
ZNDI82005	DI Water 5 Litre RFU Non-sterile	5L	Each	2 x 5L

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

RECYCLING INFORMATION

+ info see page 64



# ALPHA, BETA AND DELTA POLYCELLULOSE POUCH WIPE VERIGUARD 2

**USAGE**

Hard surface wiping

Annex 1 states that for disinfection to be effective it must be a two-stage process. The first stage should be a process to remove surface contamination and the second stage should be disinfection using a rotational system based on two types of disinfectants which have different modes of action

Micronclean has developed Delta as a sterile neutral detergent to remove surface contamination. This should be followed by Alpha and Beta as sterile rotational disinfectants which have different modes of action and prevent the long-term build-up of resistance

**BENEFITS**

- Increased liquid retention from re-sealable pouch
- Designed for small batch manufacturing
- High absorbency

**FEATURES**


- White hydro-entangled Polycellulose wipe
- 68gsm wipe weight
- 2 years shelf life
- Supplied double bagged





Product Code	Description	Size	Pack Qty	Case Qty
WSAP01009	Alpha Sessional Polycellulose C-folded Pouch Wipe Sterile	230mm x 230mm	15	10 Packs of 15
WSBT02009	Beta Sessional Polycellulose C-folded Pouch Wipe Sterile	230mm x 230mm	15	10 Packs of 15
WSND05009	Delta Sessional Polycellulose C-folded Pouch Wipe Sterile	230mm x 230mm	15	10 Packs of 15




The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



# DISPOSABLE MOP SYSTEM







RECYCLING INFORMATION

+ info see page 64



## ISOLATOR BONNET AND TOOLS

### PUREGUARD 6 AND 7

**USAGE**

Designed specifically to help clean the inside of isolators, RABs and LAF cabinets

**BENEFITS**


- Universal fit
- Dual layer polyester guarantees high levels of absorbency
- Lightweight isolator tool available in 2 handle lengths with 360° rotation

**FEATURES**

- Low linting, dual-layer 100% continuous filament polyester for increased absorbency
- 250gsm bonnet weight
- Supplied sterile one per bag ready for use
- 2 year sterile shelf life
- Handle supplied Non-Sterile ready for sterilisation

Product Code	Description	Size	Pack Qty	Case Qty
ISED0140	Dry Dual Layer Polyester Laundered Isolator Bonnet Sterile	One Size	Each	40 Packs of 1
ZNSS01035	Stainless Steel Isolator Cleaning Tool 35cm Non-Sterile	One Size	Each	1
ZNSS01061	Stainless Steel Isolator Cleaning Tool 61cm Non-Sterile	One Size	Each	1

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



## IPA AND IMS DISPOSABLE MICROFIBRE MOP STERILE PUREGUARD 1 AND 2

### USAGE

Designed for use on cleanroom floors

### BENEFITS

Uniformly dosed to achieve a coverage of 15m<sup>2</sup>  
No need for in cleanroom mixing or dosing  
2 pocket design enables a hands-free and secure connection to the mop frame  
Eliminates the risk associated with exposed chemicals and soiled water

### FEATURES

Proven 2 year shelf life  
Individually packed mops are suitable for smaller cleanrooms  
Multi packed mops are suitable for larger cleanrooms  
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
MSMI24016	IPA Individually Packed Disposable Microfibre Mop Sterile	440mm x 140mm	1	30 Packs of 1
MSMI21016	IPA Multi Packed Disposable Microfibre Mop Sterile	440mm x 140mm	10	3 Packs of 10
MSME24016	IMS Individually Packed Disposable Microfibre Mop Sterile	440mm x 140mm	1	30 Packs of 1
MSME21016	IMS Multi Packed Disposable Microfibre Mop Sterile	440mm x 140mm	10	3 Packs of 10

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



+ info see page 64



## ALPHA, BETA AND DELTA DISPOSABLE MICROFIBRE MOP STERILE PUREGUARD 4, 3 AND 5

### USAGE

Annex 1 states that for disinfection to be effective it must be a two-stage process. The first stage should be a process to remove surface contamination and the second stage should be disinfection using a rotational system based on two types of disinfectants which have different modes of action

Micronclean has developed Delta as a sterile neutral detergent to remove surface contamination. This should be followed by Alpha and Beta as sterile rotational disinfectants which have different modes of action and prevent the long-term build-up of resistance

### BENEFITS

Uniformly dosed to achieve a coverage of 15m<sup>2</sup>

No need for in cleanroom mixing or dosing

2 pocket design enables a hands-free and secure connection to the mop frame

Eliminates the risk associated with exposed chemicals and soiled water

### FEATURES

Proven 2 year shelf life

Individually packed mops are suitable for smaller cleanrooms

Multi packed mops are suitable for larger cleanrooms

Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
MSMA24016	Alpha Individually Packed Disposable Microfibre Mop Sterile	440mm x 140mm	1	30 Packs of 1
MSMA21016	Alpha Multi Packed Disposable Microfibre Mop Sterile	440mm x 140mm	10	3 Packs of 10
MSMB24016	Beta Individually Packed Disposable Microfibre Mop Sterile	440mm x 140mm	1	30 packs of 1
MSMB21016	Beta Multi Packed Disposable Microfibre Mop Sterile	440mm x 140mm	10	3 Packs of 10
MSMN24016	Delta Individually Packed Disposable Microfibre Mop Sterile	440mm x 140mm	1	30 Packs of 1
MSMN21016	Delta Multi Packed Disposable Microfibre Mop Sterile	440mm x 140mm	10	3 Packs of 10

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



RECYCLING INFORMATION



+ info see page 64

## DRY DISPOSABLE MICROFIBRE MOP STERILE

### PUREGUARD 6

#### USAGE

Designed for use on cleanroom floors

#### BENEFITS

2 pocket design enables a hands-free and secure connection to the mop frame

#### FEATURES

Proven 2 year shelf life  
Individually packed mops are suitable for smaller cleanrooms  
Multi packed mops are suitable for larger cleanrooms  
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
MSMD36016	Dry Individually Packed Disposable Microfibre Mop Sterile	440mmx 140mm	1	40 Packs of 1
MSMD36116	Dry Multi Packed Disposable Microfibre Mop Sterile	440mmx 140mm	10	4 Packs of 10

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

Please see  
individual  
product  
specifications



## RECYCLING INFORMATION



## WASTE BAGS, MOP HEADS, TROLLEY AND OTHER ACCESSORIES

### PUREGUARD 7 AND 8


#### BENEFITS

Designed to work alongside the PureGuard Mopping System



Product Code	Description	Size	Pack Qty	Case Qty
BSP1C1012	Polyethylene Bag Sterile	One Size	10	12 Packs of 10
MTSS01001	Stainless Steel Standard Mop Trolley Non-Sterile	One size	Each	1
MFSS01001	Stainless Steel Lightweight Mop Frame 40cm Non-Sterile	One Size	Each	1
MFSS79226	Stainless Steel Shrouded Locking Mop Frame 40cm Non-Sterile	One Size	Each	1
MHSS02001	Stainless Steel Fixed Mop Handle 142cm Non-Sterile	One Size	Each	1
MHSS02002	Aluminium Adjustable Mop Handle 100cm-180cm Non-Sterile	One Size	Each	1
MTSS01001	Stainless Steel Standard Mop Trolley Non-Sterile	One Size	Each	1

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

# DRY WIPES



## STERILE WIPE DECISION TREE



High Grade Cleanroom  
Medium Grade Cleanroom  
Low Grade Cleanroom

ISO Grade 4  
ISO Grade 6  
ISO Grade 8+

EU GMP Grade A  
EU GMP Grade B  
EU GMP Grade C/D





## POLYCELLULOSE WIPE STERILE WHITEGUARD 2

### USAGE

Hard surface wipe  
Transfer of product into cleanrooms

### BENEFITS

Excellent absorbency

### FEATURES

Low linting, multipurpose wipes  
68gsm wipe weight  
Hydroentangled 55% Polyester, 45% Cellulose blend

Product Code	Description	Size	Pack Qty	Case Qty
WSWG04009	Polycellulose Non-Folded Wipe Sterile	230mm x 230mm (9" x 9")	100	10 Packs of 100
WSWG06009	Polycellulose Non-Folded Wipe Sterile	230mm x 230mm (9" x 9")	300	8 Packs of 300
WSWG04012	Polycellulose Non-Folded Wipe Sterile	300mm x 300mm (12" x 12")	100	10 Packs of 100
WSCS93009	Polycellulose Non-Folded Wipe Sterile	230mm x 230mm (9" x 9")	300	8 Packs of 300
WSCS93012	Polycellulose Non-Folded Wipe Sterile	300mm x 300mm (12" x 12")	100	17 Packs of 100
WSVW11009	Polycellulose C-Folded Wipes Sterile	230mm x 230mm (9" x 9")	100	27 Packs of 100
WSVW12012	Polycellulose C-Folded Wipes Sterile	300mm x 300mm (12" x 12")	100	27 Packs of 100







RECYCLING INFORMATION

+ info see page 64



# POLYCELLULOSE WIPE NON-STERILE

## WHITEGUARD 2

**USAGE**

Hard surface wipe  
Transfer of product into cleanrooms

**BENEFITS**

Excellent absorbency

**FEATURES**

Low linting, multipurpose wipes  
68gsm wipe weight  
Hydroentangled 55% Polyester, 45% Cellulose blend

Product Code	Description	Size	Pack Qty	Case Qty
WNWG04009	Polycellulose Non-Folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	100	10 Packs of 100
WNWG06009	Polycellulose Non-Folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	300	16 Packs of 300
WNWG04012	Polycellulose Non-Folded Wipe Non-Sterile	300mm x 300mm (12" x 12")	100	20 Packs of 100
WNCS93009	Polycellulose Non-Folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	300	8 Packs of 300
WNCS93012	Polycellulose Non-Folded Wipe Non-Sterile	300mm x 300mm (12" x 12")	100	17 Packs of 100



## RECYCLING INFORMATION



+ info see page 65

## POLYESTER WIPE STERILE

### WHITEGUARD 3 AND 11

#### USAGE

Hard surface wipe  
Transfer of product into cleanrooms


#### BENEFITS


Ultra low linting wipe


#### FEATURES

Laser sealed edges to prevent linting  
100% knit continuous filament polyester  
Medium levels of absorbency

Product Code	Description	Size	Pack Qty	Case Qty
WSTG05009	Polyester Laundered Non-Folded Wipe Sterile	230mm x 230mm (9" x 9")	150	10 Packs of 150
WSDPE5010	Polyester Laundered Non-Folded Wipe Sterile	300mm x 300mm (12" x 12")	50	10 Packs of 50
WSPW80009	Polyester Laundered Non-Folded Wipe Sterile	230mm x 230mm (9"x9")	150	10 Packs of 150
MCDELN50S	Polyester Laundered Single Layer Non-Folded Wipe Sterile	300mm x 300mm (12"x12")	50	10 Packs of 50
MHPELN50S	Dual Layer Polyester Laundered Non-folded Wipe Sterile	300mm x 300mm (12" x 12")	50	4 Packs of 50







RECYCLING INFORMATION

+ info see page 65



# POLYESTER WIPE NON-STERILE

## WHITEGUARD 3

**USAGE**

Hard surface wipe  
Transfer of product into cleanrooms

**BENEFITS**

Ultra low linting wipe

**FEATURES**

Laser sealed edges to prevent linting  
100% knit continuous filament polyester  
Medium levels of absorbency

Product Code	Description	Size	Pack Qty	Case Qty
WNTG04009	Polyester Non-Folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	100	10 Packs of 100
WNTG05009	Polyester Non-Folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	150	10 Packs of 150
WNTG05012	Polyester Non-Folded Wipe Non-Sterile	300mm x 300mm (12" x 12")	150	10 Packs of 150

## RECYCLING INFORMATION



+ info see page 65



## MICRODENIER WIPE NON-STERILE WHITEGUARD 5

### USAGE

Hard surface wipe  
Excellent wipe for lens and screen cleaning


### BENEFITS

Ultra low linting wipe

### FEATURES

100% continuous filament nylon/polyester  
wipe

Product Code	Description	Size	Pack Qty	Case Qty
WNMD04009	Microdenier Non-folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	100	10 Packs of 100



RECYCLING INFORMATION

+ info see page 64



## POLYCELLULOSE ROLL WIPE NON-STERILE

### WHITEGUARD 7

#### USAGE

For use in support areas  
Hard surface wipe

#### BENEFITS

Excellent absorbency


#### FEATURES

Low linting, multipurpose wipes  
Log roll of 400 wipes

Product Code	Description	Size	Pack Qty	Case Qty
WNPW81010	Polycellulose Roll Wipe Non-Sterile White	254mm x 381mm	400 sheets/roll	2 x 1 Roll
WNPB81010	Polycellulose Roll Wipe Non-Sterile Blue	254mm x 381mm	400 sheets/roll	2 x 1 Roll

# GOGGLES

---





RECYCLING INFORMATION

+ info see page 65



## INDIVIDUALLY PACKED GOGGLE NON-STERILE VISIONGUARD 1

### USAGE

Specifically designed for use in  
Cleanroom Environments

### BENEFITS

Fully compatible with the requirements of  
Annex 1

### FEATURES


Supplied in a Non-Sterile condition but  
packed in an Autoclave bag ready for  
immediate sterilisation  
40 process lifespan for goggles without  
significant fogging  
Lightweight, supplied with adjustable head  
strap and easy-to-use snap-on clip  
Indirect ventilation to top of goggle to prevent  
fogging and protect from splashes  
Compatible with most prescription glasses  
EN 166 3 4 BT Compliant  
CE and UKCA Marked


Product Code	Description	Size	Pack Qty	Case Qty
ONACL0112	Individually Packed Goggle Non-Sterile	One Size	1	12 Packs of 1



# CLOTHING

---





RECYCLING INFORMATION

+ info see page 65



## DISPOSABLE COVERALL WITH INTEGRAL HOOD

### SUREGUARD 1

#### USAGE

Designed for cleanroom use  
Suitable for Cytotoxic work

#### BENEFITS

CE Marked to PPE Category III Type 5 & 6  
Anti-static

#### FEATURES

Single use cleanroom coverall  
CE and UKCA Marked.  
Polypropylene with polyethylene laminate fabric  
Integrated hood  
Latex free elastic at hood, wrists and ankles  
Available in Sterile and Non-Sterile Forms

Product Code	Description	Size	Pack Qty	Case Qty
DSWH04	Disposable Coverall with Integral Hood Sterile	SML - 2XL	1	25 Packs of 1
DNWH04	Disposable Coverall with Integral Hood Non-Sterile	SML - 2XL	1	25 Packs of 1



RECYCLING INFORMATION



+ info see page 65

## DISPOSABLE COVERALL WITH INTEGRAL HOOD AND FEET STERILE

### SUREGUARD 3

#### USAGE

Designed for cleanroom use  
Suitable for Cytotoxic work


#### BENEFITS

CE Marked to PPE Category III Type 5 & 6  
Anti-static

#### FEATURES

Single use cleanroom coverall  
CE and UKCA Marked.  
Polypropylene with polyethylene laminate fabric  
Integrated hood and feet  
Latex free elastic at hood, wrists and ankles

Product Code	Description	Size	Pack Qty	Case Qty
DSWH21	Disposable Coverall with Integral Hood and Feet Sterile	SML - 2XL	1	25 Packs of 1



## DISPOSABLE LABCOAT NON-STERILE SUREGUARD 2

### USAGE

Designed for laboratory and cleanroom use

### BENEFITS

CE Marked to PPE Category III Type 6  
Anti-static

### FEATURES

Single use cleanroom labcoat  
CE and UKCA Marked  
Polypropylene with polyethylene laminate fabric  
Latex free elastic at wrist

Product Code	Description	Size	Pack Qty	Case Qty
DNWH10	Disposable Labcoat Non-Sterile	SML - 2XL	1	25 Packs of 1



RECYCLING INFORMATION



+ info see page 65

## SINGLE USE CLEANROOM COMPATIBLE SOCK NON-STERILE

### TUFFGUARD 1

#### FEATURES

Specifically developed for use in cleanrooms to help with compliance with Annex 1

#### BENEFITS

Sock shaped for wearer comfort

#### FEATURES

White sock shaped  
100% synthetic fibres  
One Size

Product Code	Description	Size	Pack Qty	Case Qty
SNWH50001	Single Use Cleanroom Compatible Sock Non-Sterile	One Size	100 (50 pairs)	100 (50 pairs)



**FACEGUARD 3**  
FACEMASK WITH HEADLOOPS



**FACEGUARD 1**  
FACEMASK WITH TIES



**FACEGUARD 5 AND 6**  
BEARD SNOOD AND MOP CAP



**FACEGUARD 2**  
FACEMASK WITH EARLOOPS

## DISPOSABLE FACEMASK FACEGUARD 1, 2 AND 3 BEARD SNOOD AND MOB CAP NON-STERILE FACEGUARD 5 AND 6

### FEATURES

#### DISPOSABLE FACEMASKS

One Size

Available with Ties, Earloops or Headloops

Available in Sterile and Non-Sterile forms

### FEATURES

#### BEARD SNOOD AND MOP CAP

One Size

Comfort fit



Product Code	Description	Size	Pack Qty	Case Qty
MNTW01100	Multi Packed Tie-On Face Mask Non-Sterile	One Size	100	10 Packs of 100
MNTW01020	Individually Packed Tie-On Face Mask Non-Sterile	One Size	20	50 Packs of 20
MSTW01020	Individually Packed Tie-On Face Mask Sterile	One Size	20	15 Packs of 20
MNEW01100	Multi Packed Earloop Face Mask Non-Sterile	One Size	100	10 Packs of 100
MNEW01020	Individually Packed Earloop Face Mask Non-Sterile	One Size	20	50 Packs of 20
MSHW01020	Individually Packed Earloop Face Mask Sterile	One Size	20	25 Packs of 20
MSEW01020	Individually Packed Headloop Face Mask Sterile	One Size	20	25 Packs of 20
MNHW01100	Multi Packed Headloop Face Mask Non-Sterile	One Size	100	10 Packs of 100
MNHW01020	Individually Packed Headloop Face Mask Non-Sterile	One Size	20	50 Packs of 20
MNWH01ZZZ	Beard Snood White Non-Sterile	One Size	100	5 Packs of 100
HNWH01ZZZ	Mob Cap White Non-Sterile	One Size	100	10 Packs of 100

# GLOVES


---



## GLOVE DECISION TREE







## POLYCHLOROPRENE 300mm HAND SPECIFIC GLOVE PACKED IN L/R WALLET STERILE SKINGUARD 1

### USAGE

Specifically developed for use in cleanrooms to provide enhanced protection from cytotoxic drugs without compromising dexterity

Compatible with ISO Class 4+



### BENEFITS



Cytotoxic protection  
Enhanced Dexterity  
Latex free



### FEATURES



100% Polychloroprene  
CE / UKCA Marked  
Packed in L/R Wallet  
Hand specific  
Powder free  
Textured Palm  
Beaded cuff  
5 year shelf life

Product Code	Description	Size	Pack Qty	Case Qty
GSPG31	Polychloroprene 300mm Hand Specific Glove Packed in L/R Wallet Sterile	6.0 - 9.0	10 Pairs	5 Packs of 10 Pairs









RECYCLING INFORMATION

+ info see page 65



## LATEX 300mm HAND SPECIFIC GLOVE PACKED IN L/R WALLET STERILE SKINGUARD 5

**USAGE**

Specifically developed for use in cleanrooms to provide enhanced dexterity, flexibility and wearer comfort  
Compatible with ISO Class 4+


**BENEFITS**

Enhanced Dexterity

**FEATURES**

- 100% Latex
- CE / UKCA Marked
- Packed in L/R Wallet
- Hand specific
- Powder free
- Textured Palm
- Beaded cuff
- 5 year shelf life

Product Code	Description	Size	Pack Qty	Case Qty
GSLN05	Latex 300mm Hand Specific Glove Packed in L/R Wallet Sterile	6.0 - 9.0	10 Pairs	20 Packs of 10 Pairs



## NITRILE 300mm HAND SPECIFIC GLOVE PACKED IN L/R WALLET STERILE SKINGUARD 6

### USAGE

Specifically developed for use in cleanrooms to provide dexterity, flexibility and wearer comfort

Compatible with ISO Class 4+

### BENEFITS

Good Dexterity

Latex Free

### FEATURES

CE / UKCA Marked

Packed in L/R Wallet

Hand specific

Powder free

Textured Palm

Beaded cuff

5 year shelf life

Product Code	Description	Size	Pack Qty	Case Qty
GSNW06	Nitrile 300mm Hand Specific Glove Packed in L/R Wallet Sterile	6.0 - 9.0	10 Pairs	20 Packs of 10 Pairs



RECYCLING INFORMATION

+ info see page 65



## NITRILE 405mm HAND SPECIFIC GLOVE PACKED IN L/R WALLET STERILE SKINGUARD 7

### USAGE

Specifically developed for use in cleanrooms to provide dexterity, flexibility and wearer comfort

Compatible with ISO Class 4+

### BENEFITS

Good Dexterity

Latex Free

### FEATURES

CE / UKCA Marked

Packed in L/R Wallet

Hand specific

Powder free

Textured Palm

Beaded cuff

5 year shelf life

Product Code	Description	Size	Pack Qty	Case Qty
GSNW07	Nitrile 405mm Hand Specific Glove Packed in L/R Wallet Sterile	6.0 - 9.0	10 Pairs	20 Packs of 10 Pairs



RECYCLING INFORMATION



+ info see page 65

## NITRILE 230mm AMBIDEXTROUS GLOVE BLUE BOXED NON-STERILE SKINGUARD 2

### USAGE

Specifically developed for use in cleanrooms  
Compatible with ISO Class 7/8

### BENEFITS

Latex Free  
Bulk packed in cardboard dispensers

### FEATURES

Blue  
Powder free  
Nitrile material  
Ambidextrous  
Textured surface for increased grip  
Beaded cuff

Product Code	Description	Size	Pack Qty	Case Qty
GNNB02	Nitrile 230mm Ambidextrous Glove Blue Boxed Non-Sterile	SML - XLG	150 Gloves	10 Packs of 150



+ info see page 65



## LATEX 300mm AMBIDEXTROUS GLOVE BULK PACKED NON-STERILE SKINGUARD 10

### USAGE

Specifically developed for use in cleanrooms to provide enhanced dexterity, flexibility and wearer comfort  
Compatible with ISO Class 4+

### BENEFITS

Enhanced Dexterity


### FEATURES

100% Latex  
CE / UKCA Marked  
Powder free  
Textured Palm  
Beaded cuff  
Bulk Packed

Product Code	Description	Size	Pack Qty	Case Qty
GNLN10	Latex 300mm Ambidextrous Glove Bulk Packed Non-Sterile	SML - XLG	100 Gloves	10 Packs of 100



RECYCLING INFORMATION



+ info see page 65

## NITRILE 300mm AMBIDEXTROUS GLOVE BULK PACKED NON-STERILE SKINGUARD 11

### USAGE

Specifically developed for use in cleanrooms to provide dexterity, flexibility and wearer comfort  
Compatible with ISO Class 4+

### BENEFITS

Good Dexterity  
Latex Free  
Liquid particle count <600 at 0.5µm

### FEATURES

CE / UKCA Marked  
Powder free  
Textured Palm  
Beaded cuff  
Bulk Packed

Product Code	Description	Size	Pack Qty	Case Qty
GNNW11	Nitrile 300mm Ambidextrous Glove Bulk Packed Non-Sterile	SML - XLG	100 Gloves	10 Packs of 100



+ info see page 65



## NITRILE 300mm AMBIDEXTROUS GLOVE BULK PACKED NON-STERILE SKINGUARD 12

### USAGE

Specifically developed for use in cleanrooms to provide dexterity, flexibility and wearer comfort  
Compatible with ISO Class 6+

### BENEFITS

Good Dexterity  
Latex Free  
Liquid particle count <3000 at 0.5µm


### FEATURES


CE / UKCA Marked  
Powder free  
Textured Palm  
Beaded cuff  
Bulk Packed

Product Code	Description	Size	Pack Qty	Case Qty
GNNN12	Nitrile 300mm Ambidextrous Glove Bulk Packed Non-Sterile	SML- XLG	100 Gloves	10 Packs of 100



# STATIONERY





**GUARDMARK 1**  
CLEANROOM BALLPOINT PEN STERILE



**GUARDMARK 2**  
CLEANROOM MARKER PEN BLACK STERILE

## CLEANROOM BALLPOINT PEN STERILE GUARDMARK 1

## CLEANROOM MARKER PEN BLACK STERILE GUARDMARK 2

### USAGE

#### CLEANROOM BALLPOINT PEN

Designed for controlled environments, this sterile ballpoint pen features low smear, low sodium black ink and inherent ESD properties

#### CLEANROOM MARKER PEN

Proven ability to write on most typical laboratory surfaces including plastic, glass, and stainless steel

### BENEFITS

#### CLEANROOM BALLPOINT PEN

Low smear black ink

#### CLEANROOM MARKER PEN

Low particulate

### FEATURES

#### CLEANROOM BALLPOINT PEN

Low smear black ink with low sodium content  
Inherent ESD properties  
Packed in ISO 14644-1 ISO Class 6 cleanroom  
2 years shelf life

#### CLEANROOM MARKER PEN

Quick-drying ink which is water and UV-light resistant  
Does not dry out even if the cap is left off for a day  
Packed in ISO 14644-1 ISO Class 6 cleanroom  
2 years shelf life

Product Code	Description	Size	Pack Qty	Case Qty
ZSBK11210	Cleanroom Ballpoint Pen Sterile	One Size	12	10 Packs of 12
ZSBK03010	Cleanroom Marker Pen Black Sterile	One Size	1	10 Packs of 1



## RECYCLING INFORMATION



+ info see page 66

## CLEANROOM PAPER NON-STERILE

### GUARDMARK 3

#### USAGE

Designed for in cleanroom use with laser or inkjet printers

#### BENEFITS

Available in 5 colours: White, Blue, Pink, Green and Yellow

Low particulate generation

#### FEATURES

250 sheets per pack

Basic weight 72gsm

Product Code	Description	Size	Pack Qty	Case Qty
ZNWH04072	Cleanroom Paper White Non-Sterile	A4 - 18lbs	250 sheets	10 Packs of 250 sheets
ZNBL04072	Cleanroom Paper Blue Non-Sterile	A4 - 18lbs	250 sheets	10 Packs of 250 sheets
ZNPI04072	Cleanroom Paper Pink Non-Sterile	A4 - 18lbs	250 sheets	10 Packs of 250 sheets
ZNGR04072	Cleanroom Paper Green Non-Sterile	A4 - 18lbs	250 sheets	10 Packs of 250 sheets
ZNYE04072	Cleanroom Paper Yellow Non-Sterile	A4 - 18lbs	250 sheets	10 Packs of 250 sheets



## CLEANROOM MUNISING PAPER STERILE AND NON-STERILE GUARDMARK 4

### USAGE

Designed for in cleanroom use with laser or inkjet printers

### BENEFITS

Low particulate generation  
Good heat stability  
ROHS Compliant


### FEATURES


Available in white, yellow and blue  
Manufactured from a combination of Cellulose Fibre and Synthetic Latex  
250 sheets per pack  
Basic weight 80gsm

Product Code	Description	Size	Pack Qty	Case Qty
ZNWH04022	Cleanroom Munising Paper White Non-Sterile	A4	250 sheets	10 Packs of 250 sheets
ZNBL04022	Cleanroom Munising Paper Blue Non-Sterile	A4	250 sheets	10 Packs of 250 sheets
ZNYE04022	Cleanroom Munising Paper Yellow Non-Sterile	A4	250 sheets	10 Packs of 250 sheets
ZSWH04225	Cleanroom Munising Paper White Sterile	A4	250 sheets	10 Packs of 250 sheets
ZSGR04022	Cleanroom Munising Paper Green Sterile	A4	250 sheets	10 Packs of 250 sheets
ZSBL04225	Cleanroom Munising Paper Blue Sterile	A4	250 sheets	10 Packs of 250 sheets

# OVERSHOES

---





## POLYTHENE OVERSHOE BLUE NON-STERILE

### TUFFGUARD 2

#### FEATURES

Elastic opening for easy donning

Available in two sizes

Available in two materials depending on usage

Manufactured from PE

Product Code	Description	Size	Pack Qty	Case Qty
SNBL01016	Blue Heavy Duty Polythene Overshoe 16" Non-Sterile	420mm x 160mm	100	10 Packs of 100
SNBL31014	Blue Polythene Overshoe 14" Non-Sterile	360mm x 150mm	100	20 Packs of 100
SNBL31016	Blue Polythene Overshoe 16" Non-Sterile	410mm x 150mm	100	20 Packs of 100



RECYCLING INFORMATION



+ info see page 66

## PVC OVERSHOE NON-STERILE

### TUFFGUARD 3

#### FEATURES

- Non slip sole
- Available in two widths and two lengths
- Manufactured from PVC
- Elastic opening for easy donning

Product Code	Description	Size	Pack Qty	Case Qty
SNWH81006	PVC 6"x14" Overshoe White Non-Sterile	150 mm × 355 mm	500	500
SNWH83016	PVC 6"x16" Overshoe White Non-Sterile	150 mm × 405 mm	500	500
SNWH82009	PVC 9"x14" Overshoe White Non-Sterile	230 mm × 355 mm	500	500
SNWH81016	PVC 9"x16" Overshoe White Non-Sterile	230 mm × 405 mm	500	500
SNRE81009	PVC 9"x16" Overshoe Red Non-Sterile	230 mm × 405 mm	500	500



+ info see page 66



## OVERSHOE WITH TEXTURED PVC SOLE NON-STERILE TUFFGUARD 4

## OVERSHOE ANTI-SLIP NON-STERILE TUFFGUARD 6

### BENEFITS

#### TUFFGUARD 4

Stitched seams with sole

#### TUFFGUARD 6

Anti slip sole

### FEATURES

#### TUFFGUARD 4

Microporous polypropylene with polyethylene laminate fabric  
Polyvinylchloride (PVC) soles


#### TUFFGUARD 6

Elastic opening for easy donning  
Manufactured from Polypropylene (+EVA) laminate

Product Code	Description	Size	Pack Qty	Case Qty
SNWH87016	Overshoe with textured PVC sole Non-sterile	180mm x 400mm	100	10 Packs of 100
SNWH03016	Overshoe White Anti-Slip Non-Sterile	125mm x 420mm	100	10 Packs of 100



# MATS





+ info see page 66



## ACCESS MAT MEDIUM TACK NON-STERILE GUARDTACK 1

### USAGE

Peel-Off access mat designed to capture contamination before entry into a controlled area

### BENEFITS

Medium "Tack" formulation gives a higher level of dirt capture than GuardTack 2 whilst not inhibiting free movement

### FEATURES

30 Individually numbered LDPE peel-off sheets per contamination control mat  
8 mats per pack  
Available in three colours (Blue / Grey / White)  
Available in three sizes

Product Code	Description	Size	Pack Qty	Case Qty
ANXX08018	Contamination Control Mats Medium Tack Non-Sterile	458mm x 1145mm (18" x 45")	30 sheets/mat	8 Mats
ANXX08026	Contamination Control Mats Medium Tack Non-Sterile	660mm x 1145mm (26" x 45")	30 sheets/mat	8 Mats
ANXX08336	Contamination Control Mats Medium Tack Non-Sterile	914mm x 914mm (36" x 36")	30 sheets/mat	8 Mats

Where XX = BL for Blue, GY for Grey and WH for White



## ACCESS MAT LOW TACK NON-STERILE GUARDTACK 2

### USAGE

Peel-Off access mat designed to capture contamination before entry into a controlled area

### BENEFITS

Low "Tack" formulation gives a lower level of dirt capture compared to GuardTack 1 but does not affect free movement

### FEATURES


30 Individually numbered LDPE peel-off sheets per contamination control mat  
4 mats per pack  
Available in three colours (Blue / Grey / White)  
Available in five sizes

Product Code	Description	Size	Pack Qty	Case Qty
ANXX04018	Contamination Control Mats Low Tack Non-Sterile	458mm x 1145mm (18" x 45")	30 sheets/mat	4 Mats
ANXX04236	Contamination Control Mats Low Tack Non-Sterile	660mm x 914mm (24" x 36")	30 sheets/mat	4 Mats
ANXX04026	Contamination Control Mats Low Tack Non-Sterile	660mm x 1145mm (26" x 45")	30 sheets/mat	4 Mats
ANXX04036	Contamination Control Mats Low Tack Non-Sterile	914mm x 1145mm (36" x 45")	30 sheets/mat	4 Mats
ANXX04336	Contamination Control Mats Low Tack Non-Sterile	914mm x 914mm (36" x 36")	30 sheets/mat	4 Mats

Where XX = BL for Blue, GY for Grey and WH for White



+ info see page 66



## BLUE WALK ACCESS MAT WASHABLE NON-STERILE GUARDTACK 3

### USAGE

Washable, high performance contamination control mat for use with medium traffic loads


### BENEFITS

Specially designed and manufactured to capture, retain and reduce foot and wheel-borne viable and non-viable particulates, over an extended period of time

### FEATURES

Installed by trained service technicians.  
3-5 year in service life  
Silver based antimicrobial embedded into the surface  
Max 1000kg per mat area

Product Code	Description	Size	Pack Qty	Case Qty
ANBLWLK12	Blue Walk Access Mat 1.2m Non-Sterile	1200mm	Each	1 Mat
ANBLWLK20	Blue Walk Access Mat 2.0m Non-Sterile	2000mm	Each	1 Mat



## GREY PRO ACCESS MAT WASHABLE NON-STERILE

### GUARDTACK 4

#### USAGE

Washable, high performance contamination control mat for use with high traffic loads

#### BENEFITS

Specially designed and manufactured to capture, retain and reduce foot and wheel-borne viable and non-viable particulates, over an extended period of time


#### FEATURES

Installed by trained service technicians  
 3-5 year in service life  
 Silver based antimicrobial embedded into the surface  
 Max 3000kg per mat area  
 Available with stainless steel edging to help protect the entry and exit off the mat surface

Product Code	Description	Size	Pack Qty	Case Qty
ANGYPRO12	Grey Pro Access Mat 1.2m Non-Sterile	1200mm	Each	1 Mat
ANGYPRO20	Grey Pro Access Mat 2.0m Non-Sterile	2000mm	Each	1 Mat



+ info see page 66



## BLUE SELF-INSTALL ACCESS MAT NON-STERILE GUARDTACK 5

### USAGE

Washable, high performance contamination control mat for use with light traffic loads

### BENEFITS

Specially designed and manufactured to capture, retain and reduce foot and wheel-borne viable and non-viable particulates, over an extended period of time

### FEATURES

Available in 3 sizes  
Self-install  
3-5 year in service life  
Silver based antimicrobial embedded into the surface  
Max 400kg per mat area

Product Code	Description	Size	Pack Qty	Case Qty
ANBLINS06	Blue Self-install Access Mat 0.6m x 1.2m Non-Sterile	1200mm x 600mm	Each	1 Mat
ANBLINS10	Blue Self-install Access Mat 1.0m x 1.2m Non-Sterile	1000mm x 1200mm	Each	1 Mat
ANBLINS15	Blue Self-install Access Mat 1.0m x 1.5m Non-Sterile	1000mm x 1500mm	Each	1 Mat

# MICRONGREEN





**“Our aim is to be the most sustainable cleanroom textile and cleanroom consumables business on the planet.”**


Micronclean is committed to contributing to a sustainable future based on the three pillars of:

- Protecting our climate
- Creating our circular economy
- Delivering a positive impact for our people and communities





# RECYCLING INFORMATION



## RECYCLING INFORMATION



### HAND SANITISER

Page 01

Component Name	Component Material	Disposal
Bottle Pump	Mixed Materials	Non-recyclable
500ml Bottle	PET	Recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### TRIGGER SPRAY

Page 05 | 12

Component Name	Component Material	Disposal
Trigger Spray Head with Dip Tube	Mixed Materials	Non-recyclable
500ml/950ml Bottle	HDPE	Recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### 5 LITRE RFU

Page 06 | 13

Component Name	Component Material	Disposal
Container Cap	HDPE	Recyclable
5L Container	PP	Recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### POLYPROPYLENE TUB WIPES

Page 07

Component Name	Component Material	Disposal
Wipe	Polypropylene	Recyclable
Tub	HDPE	Recyclable
Tub Lid	PE	Recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### POLYCELLULOSE POUCH WIPES

Page 08 | 14

Component Name	Component Material	Disposal
Wipe	Polycellulose	Non-recyclable
Wipe Pouch	Polyester and Polyethylene blend	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### IPA POLYESTER POUCH WIPES

Page 09

Component Name	Component Material	Disposal
Wipe	Polyester	Recyclable
Wipe Pouch	Polyester and Polyethylene blend	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### IMS SILVER PACK POLYETHYLENE/POLYPROPYLENE WIPES

Page 10

Component Name	Component Material	Disposal
Wipe	Polyethylene / Polypropylene	Non-recyclable
Wipe Pouch	OPP + PE + VMPET	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### ISOLATOR BONNET AND TOOLS

Page 16

Component Name	Component Material	Disposal
Bonnet	Polyester	Recyclable
Handle	Stainless Steel	Recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### DISPOSABLE MICROFIBRE MOP STERILE

Page 17 | 18 | 19

Component Name	Component Material	Disposal
Mop	Microfibre (Polyester and Polyamide)	Non-recyclable
Mop Pouch (Individually Packed Only)	Polyester and Polyethylene blend	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### POLYCELLULOSE WIPES / ROLL WIPES

Page 23 | 24 | 28

Component Name	Component Material	Disposal
Wipe	Hydroentangled 55% Polyester, 45% Cellulose blend	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable

## RECYCLING INFORMATION



### POLYESTER WIPES STERILE

Page 25 | 26

Component Name	Component Material	Disposal
Wipe	100% Polyester	Recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### MICRODENIER WIPES NON-STERILE

Page 27

Component Name	Component Material	Disposal
Wipe	100% continuous filament nylon/ polyester wipe	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### INDIVIDUALLY PACKED GOGGLE NON-STERILE

Page 30

Component Name	Component Material	Disposal
Goggle	Thermoplastic Rubber, Silicone, Polycarbonate	Non-recyclable
Autoclave Bag	French Paper With Blue laminated Film (PE and CPP blend)	Non-recyclable
Outer Bag	Nylon	Recyclable
Box	Cardboard	Recyclable



### DISPOSABLE COVERALL / DISPOSABLE LAB COAT

Page 32 | 33 | 34

Component Name	Component Material	Disposal
Coverall	White Microporous PP with PE laminated with zip and rubber elastic	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### CLEANROOM POLYESTER SOCKS

Page 35

Component Name	Component Material	Disposal
Socks	Mixed synthetic fibre	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### DISPOSABLE FACEMASKS, BEARD SNOOD AND MOB CAP

Page 36

Component Name	Component Material	Disposal
Mobcap/Beard Snood	PP Fabric With Rubber Elastic	Non-recyclable
Facemask	Combined polymer ES (Ethylene Side By Side), Melt Blown PP (MBPP), Polyurethane, metal stripped with coating	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### POLYCHLOROPRENE GLOVES

Page 39

Component Name	Component Material	Disposal
Gloves	Polychloroprene	Recyclable
L & R Wallet	HDPE	Recyclable
Polybag	HDPE	Recyclable
Box	Cardboard	Recyclable



### LATEX GLOVES

Page 40 | 44

Component Name	Component Material	Disposal
Gloves	Latex	Recyclable
L & R Wallet	HDPE	Recyclable
Polybag	MDPE	Recyclable
Box	Cardboard	Recyclable



### NITRILE GLOVES

Page 41 | 42 | 43 | 45 | 46

Component Name	Component Material	Disposal
Gloves	Nitrile Butadiene Rubber	Recyclable
L & R Wallet	HDPE	Recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### CLEANROOM BALLPOINT PEN AND MARKER PEN BLACK

Page 48

Component Name	Component Material	Disposal
Pen	Polypropylene / Polyester	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable

## RECYCLING INFORMATION



### CLEANROOM PAPER

Page 49

Component Name	Component Material	Disposal
Paper	Paper	Recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### PVC OVERSHOE / PVC SOLE / OVERSHOE ANTI-SLIP

Page 53 | 54

Component Name	Component Material	Disposal
Overshoe	PVC With Rubber Elastic Band	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### CLEANROOM MUNISING PAPER

Page 50

Component Name	Component Material	Disposal
Paper	Cellulose Fiber and Synthetic Latex	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### ACCESS MAT MEDIUM AND LOW TACK

Page 56 | 57

Component Name	Component Material	Disposal
Mat	LDPE coated with water-based acrylic adhesive	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### POLYTHENE OVERSHOE BLUE

Page 52

Component Name	Component Material	Disposal
Overshoe (Heavy Duty)	Chlorinated PE (CPE) With A Rubber Elastic Band	Non-recyclable
Overshoe	LDPE/ LLDPE with a Rubber Elastic Band	Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable



### ACCESS MAT


Page 58 | 59 | 60

Component Name	Component Material	Disposal
Mat		Non-recyclable
Polybag	LDPE	Recyclable
Box	Cardboard	Recyclable

#### Recycling information relates to

- UK Law and Practices
- Clean items that are free from contamination and residues

# MICRONCLEAN ASSEMBLY



## MICRONCLEAN ASSEMBLY



The current Micronclean manufacturing cleanroom for cleanroom consumables was commissioned in 2011 and has a dedicated manufacturing area of 61m<sup>2</sup>. When the facility was built this was perfect for the cleanroom consumables market we saw at that time and with various improvement to the equipment we have improved the production capacity to 1,350 boxes a week on average.

However, since we commissioned the unit, we have seen sustained growth in demand and since 2020 we have seen double digit growth year on year which means we have outgrown the original production facility even with its enhancements. As a result, in November of 2020 we committed £1.8million to build and equip a brand-new cleanroom manufacturing facility, Micronclean Assembly with a dedicated manufacturing area of 284m<sup>2</sup>.

From the outset Micronclean Assembly was designed to be a paper free environment and as a result has become our second paperless manufacturing facility.

Inside the cleanroom the quadrupling of manufacturing area has allowed Micronclean to install three dedicated and distinct process lines:



**Wet Line 1** is a fully automated filling and sealing line dedicated to the production of trigger sprays and tub wipes. The line has the capacity to process 20 trigger spray bottles or tubs a minute.



**Wet Line 2** is a semi-automatic filling and sealing line dedicated to the production of pouch wipes and disposable mops.



**Dry Line 1** is a manual packing line for dry good with the facility to vacuum pack as required so is ideal for isolator bonnets, dry wipes and coveralls.

The installation of the three new lines has unlocked a large amount of potential which can be detailed as follows:

- A large increase in capacity.
- Improved quality control as we move away from very manual processes to either automatic or semi-automatic processing.

This all leads to increased process capability and as a result as a Micronclean Assembly is delivering greater reliability of supply and timing.

---

# NOTES

---



Roman Bank | Skegness | Lincolnshire | PE25 1SQ | UK  
T: +44 (0)1754 767 377  
E: [sales@micronclean.co.uk](mailto:sales@micronclean.co.uk)  
W: [www.micronclean.com](http://www.micronclean.com)



MCB 05/2025