



PRODUCT SPECIFICATION

REUSABLE CLEANROOM SOCK

PROTRACK CODE | SO5

PRODUCT INFORMATION

PRODUCT & SERVICE FEATURES

The Micronclean Reusable Sock service is tailored to your needs and is;

- Designed specifically for the cleanroom environment.
- Fully compliant with the requirements detailed in EU GMP Annex 1.
- Able to offer a major reduction in waste going to landfill compared with single use disposable socks.
- Fully RFID tagged so that it is compatible with our Protrack garment tracking system.

Alongside this the socks itself has not been forgotten in that it offers socks that;

- Are logo'd to support simple visual checks that show if the correct socks are being worn.
- Generate increased wearer comfort by being 'sock' shaped.
- Are available in three colour coded sizes to again generate long term wearer comfort.
- Feature a COOLMAX® Technical Fabric giving;
 - A soft feel
 - High Levels of Breathability
 - Superb Moisture Wicking

PRODUCT DESCRIPTION

Micronclean offer a Laundered Reuseable Cleanroom Sock service, allowing our ISO Class 4 customers amongst others to be fully compliant with EU GMP Annex 1.

MANUFACTURING

Processed and packaged in an ISO Class B (GMP Grade B) Cleanroom Laundry in accordance in RABC principles.

PACKAGING

Socks are packed into single sized sealed polythene packs of up to 50 socks.

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems.

STORAGE CONDITIONS

Store in cool, well ventilated area.

MATERIAL DIMENSIONS

Sock Material:

- 65% Synthetic Technical Fibre
- 22% Polyester
- 3% Rubber
- 6% Nylon
- 4% Elastic Fibre

Sock Size:

- Small - Yellow Band (up to size 6)
- Medium - Red Band (size 6 to 10)
- Large - Blue Band (size 10 to 13)