



## PRODUCT SPECIFICATION

### STERILE PACKS

## BD 10ML

### SYRINGE MULTIPACKS

#### PRODUCT INFORMATION

##### PRODUCT DESCRIPTION

Three piece BD single use syringe packs

##### PRODUCT FEATURES

###### Materials

Barrel and plunger – polypropylene  
Gasket – latex-free, medical grade rubber  
Lubricant – medical grade silicone oil  
( $<0.25 \text{ mg/cm}^2$ )

###### Connection type

Luer-Lock (central)

###### Applicable standards

ISO 7886-1 / ISO 7886-2

###### Medical device classification

Class I (sterile)

##### MANUFACTURING

Becton Dickinson & Company 1 Becton  
Drive Franklin Lakes NJ 07417 USA

##### PACKAGING

Four layer laminate  
Easy tear notch  
Alcohol and VHP compatible

##### STERILISATION DETAILS

Gamma irradiation min. 25kGy  
Complies with applicable parts of ISO 11137  
series

##### QUALITY ASSURANCE

Manufactured under a BS EN ISO  
13485:2016 quality management system

##### SHELF LIFE

3 years from date of manufacture

##### STORAGE CONDITIONS

Store in a dry, well ventilated area and  
protect from direct sunlight and heat  
sources

## TECHNICAL INFORMATION

Compatible with gamma irradiation  
 BSE/TSE – considered not to present any risk\*  
 Phthalate free  
 Latex free  
 PVC free  
 Single scale easy to read graduation

\*raw materials may contain tallow-derived substances which have been produced under conditions which meet or exceed the requirements of EN 22442-1 and EMA 410/01 and are therefore considered not to present any risk in respect of TSE/BSE or other animal-borne diseases.

## ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
10ML05CC	10ml BD Syringe Sterile 10 Packs of 5	10ml	5	50	280 x 170 x 200mm	0.70 KG
10ML10CC	10ml BD Syringe Sterile 10 Packs of 10	10ml	10	100	280 x 170 x 390mm	1.30 KG
10ML25CC	10ml BD Syringe Sterile 10 Packs of 25	10ml	25	250	280 x 340 x 400mm	2.90 KG
10ML25CH	10ml BD Syringe and Cap Sterile 10 Packs of 25	10ml	25	250	280 x 340 x 400 mm	2.90 KG
10ML25SC	10ml BD Syringe and Cap Sleeved Sterile 8 Packs of 25	10ml	25	200	560 x 330 x 380 mm	3.70 KG
10ML50CC	10ml BD Syringe Sterile 5 Packs of 50	10ml	50	250	280 x 340 x 400 mm	2.90 KG