



## PRODUCT SPECIFICATION STERILE PACKS

## BD 1ML SYRINGE MULTIPACKS

### PRODUCT INFORMATION

#### PRODUCT DESCRIPTION

Three piece BD single use syringe packs

#### PRODUCT FEATURES

##### Materials

Barrel and plunger – polycarbonate  
Gasket – latex-free, medical grade rubber  
Lubricant – medical grade silicone oil  
( $<0.25 \text{ mg/cm}^2$ )

##### Connection type

Luer-Lock (central)

##### Applicable standards

ISO 7886-1 / ISO 7886-2  
ISO 80369-7

##### Medical device classification

Class I (sterile)

#### MANUFACTURING

Becton Dickinson & Company 1 Becton  
Drive Franklin Lakes NJ 07417 USA

#### PACKAGING

Four layer laminate  
Easy tear notch  
Alcohol and VHP compatible

#### STERILISATION DETAILS

Gamma irradiation min. 25kGy  
Complies with applicable parts of ISO 11137  
series

#### QUALITY ASSURANCE

Manufactured under a BS EN ISO  
13485:2016 quality management system

#### SHELF LIFE

3 years from date of manufacture

#### STORAGE CONDITIONS

Store in a dry, well ventilated area and  
protect from direct sunlight and heat  
sources

## TECHNICAL INFORMATION

Compatible with gamma irradiation  
BSE/TSE – considered not to present any risk\*  
Phthalate free  
Latex free  
PVC free  
Contains Bisphenol A  
Single scale easy to read graduation

\*raw materials may contain tallow-derived substances which have been produced under conditions which meet or exceed the requirements of EN 22442-1 and EMA 410/01 and are therefore considered not to present any risk in respect of TSE/BSE or other animal-borne diseases.

## ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
01ML05CC	1ml luer lock syringe plus tip protector, 5 pack	1ml	5	100	280 x 170 x 200mm	1.00 KG
01ML10CC	1ml luer lock syringe plus tip protector, 10 pack	1ml	10	100	280 x 170 x 200mm	1.00 KG
01ML25CC	1ml luer lock syringe plus tip protector, 25 pack	1ml	25	250	280 x 170 x 200mm	2.00 KG
01ML25SC	1ml luer lock syringe plus barcode and caps, 12 x 25 pack	1ml	25	300	280 x 340 x 400 mm	4.30 KG