



## PRODUCT SPECIFICATION

### STERILE PACKS

## BD 20ML

### SYRINGE MULTIPACKS

#### PRODUCT INFORMATION

##### PRODUCT DESCRIPTION

Three piece BD single use syringe packs

##### PRODUCT FEATURES

###### Materials

Barrel and plunger – polypropylene

Gasket – latex-free, medical grade rubber

Lubricant – medical grade silicone oil  
( $<0.25 \text{ mg/cm}^2$ )

###### Connection type

Luer-Lock (central)

###### Applicable standards

ISO 7886-1 / ISO 7886-2

ISO 80369-7

###### Medical device classification

Class IIa

##### MANUFACTURING

Becton Dickinson & Company 1 Becton  
Drive Franklin Lakes NJ 07417 USA

##### PACKAGING

Four layer laminate

Easy tear notch

Alcohol and VHP compatible

##### STERILISATION DETAILS

Gamma irradiation min. 25kGy

Complies with applicable parts of ISO 11137  
series

##### QUALITY ASSURANCE

Manufactured under a BS EN ISO  
13485:2016 quality management system

##### SHELF LIFE

3 years from date of manufacture

##### STORAGE CONDITIONS

Store in a dry, well ventilated area and  
protect from direct sunlight and heat  
sources

## TECHNICAL INFORMATION

Compatible with gamma irradiation  
BSE/TSE – considered not to present any risk\*  
Phthalate free  
Latex free  
PVC free  
Single scale easy to read graduation

\*raw materials may contain tallow-derived substances which have been produced under conditions which meet or exceed the requirements of EN 22442-1 and EMA 410/01 and are therefore considered not to present any risk in respect of TSE/BSE or other animal-borne diseases.

## ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
20ML05CC	20ml luer lock syringe plus tip protector, 5 pack	20ml	5	50	280 x 170 x 390mm	0.90 KG
20ML10CC	20ml luer lock syringe plus tip protector, 10 pack	20ml	10	100	280 x 170 x 390mm	1.80 KG
20ML25CC	20ml luer lock syringe plus tip protector, 25 pack	20ml	25	250	280 x 340 x 400mm	4.40 KG
20ML50CC	20ml luer lock syringe plus tip protector, 50 pack	20ml	50	500	560 x 330 x 380 mm	8.50 KG
20ML25SC	20ml luer lock syringe plus barcode and cap, sleeved, 8 x 25 pack	20ml	25	200	560 x 330 x 380 mm	5.50 KG