



PRODUCT SPECIFICATION

STERILE PACKS

BD 20ML

SYRINGE MULTIPACKS

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Three piece BD single use syringe packs

PRODUCT FEATURES

Materials

Barrel and plunger – polypropylene
Gasket – latex-free, medical grade rubber
Lubricant – medical grade silicone oil
($<0.25 \text{ mg/cm}^2$)

Connection type

Luer-Lock (central)

Applicable standards

ISO 7886-1 / ISO 7886-2

Medical device classification

Class IIa

MANUFACTURING

Becton Dickinson & Company 1 Becton
Drive Franklin Lakes NJ 07417 USA

PACKAGING

Four layer laminate
Easy tear notch
Alcohol and VHP compatible

STERILISATION DETAILS

Gamma irradiation min. 25kGy
Complies with applicable parts of ISO 11137
series

QUALITY ASSURANCE

Manufactured under a BS EN ISO
13485:2016 quality management system

SHELF LIFE

3 years from date of manufacture

STORAGE CONDITIONS

Store in a dry, well ventilated area and
protect from direct sunlight and heat
sources

TECHNICAL INFORMATION

Compatible with gamma irradiation
 BSE/TSE – considered not to present any risk*
 Phthalate free
 Latex free
 PVC free
 Single scale easy to read graduation

*raw materials may contain tallow-derived substances which have been produced under conditions which meet or exceed the requirements of EN 22442-1 and EMA 410/01 and are therefore considered not to present any risk in respect of TSE/BSE or other animal-borne diseases.

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
20ML05CC	20ml BD Syringe Sterile 10 Packs of 5	20ml	5	50	280 x 170 x 390mm	0.90 KG
20ML10CC	20ml BD Syringe Sterile 10 Packs of 10	20ml	10	100	280 x 170 x 390mm	1.80 KG
20ML25CC	20ml BD Syringe Sterile 10 Packs of 25	20ml	25	250	280 x 340 x 400mm	4.40 KG
20ML50CC	20ml BD Syringe Sterile 10 Packs of 50	20ml	50	500	560 x 330 x 380 mm	8.50 KG
20ML25SC	20ml BD Syringe and Cap Sleeved Sterile 8 Packs of 25	20ml	25	200	560 x 330 x 380 mm	5.50 KG