



PRODUCT SPECIFICATION

STERILE PACKS

BD 30ML

SYRINGE MULTIPACKS

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Three piece BD single use syringe packs

PRODUCT FEATURES

Materials

Barrel and plunger – polypropylene

Gasket – latex-free, medical grade rubber

Lubricant – medical grade silicone oil
($<0.25 \text{ mg/cm}^2$)

Connection type

Luer-Lock (central)

Applicable standards

ISO 7886-1 / ISO 7886-2

ISO 80369-7

Medical device classification

Class IIa

MANUFACTURING

Becton Dickinson & Company 1 Becton
Drive Franklin Lakes NJ 07417 USA

PACKAGING

Four layer laminate

Easy tear notch

Alcohol and VHP compatible

STERILISATION DETAILS

Gamma irradiation min. 25kGy

Complies with applicable parts of ISO 11137
series

QUALITY ASSURANCE

Manufactured under a BS EN ISO
13485:2016 quality management system

SHELF LIFE

3 years from date of manufacture

STORAGE CONDITIONS

Store in a dry, well ventilated area and
protect from direct sunlight and heat
sources

TECHNICAL INFORMATION

Compatible with gamma irradiation
BSE/TSE – considered not to present any risk*
Phthalate free
Latex free
PVC free
Single scale easy to read graduation

*raw materials may contain tallow-derived substances which have been produced under conditions which meet or exceed the requirements of EN 22442-1 and EMA 410/01 and are therefore considered not to present any risk in respect of TSE/BSE or other animal-borne diseases.

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
30ML05CC	30ml luer lock syringe plus tip protector, 5 pack	30ml	5	50	280 x 340 x 400 mm	1.70 KG
30ML10CC	30ml luer lock syringe plus tip protector, 10 pack	30ml	10	100	280 x 340 x 400 mm	2.70 KG
30ML25CC	30ml luer lock syringe plus tip protector, 25 pack	30ml	25	250	560 x 330 x 380 mm	6.00 KG
30ML25CH	30ml luer lock syringe plus tip protector and cap, 50 pack	30ml	25	250	560 x 330 x 380 mm	6.00 KG
30ML25SC	30ml luer lock syringe plus barcode and cap, sleeved, 8 x 25 pack	30ml	25	200	560 x 330 x 380 mm	6.40 KG