



## PRODUCT SPECIFICATION

### STERILE PACKS

## BD 3ML

### SYRINGE MULTIPACKS

#### PRODUCT INFORMATION

##### PRODUCT DESCRIPTION

Three piece BD single use syringe packs

##### PRODUCT FEATURES

###### Materials

Barrel and plunger – polypropylene

Gasket – latex-free, medical grade rubber

Lubricant – medical grade silicone oil  
( $<0.25 \text{ mg/cm}^2$ )

###### Connection type

Luer-Lock (central)

###### Applicable standards

ISO 7886-1 / ISO 7886-2

ISO 80369-7

###### Medical device classification

Class I (sterile)

##### MANUFACTURING

Becton Dickinson & Company 1 Becton  
Drive Franklin Lakes NJ 07417 USA

##### PACKAGING

Four layer laminate

Easy tear notch

Alcohol and VHP compatible

##### STERILISATION DETAILS

Gamma irradiation min. 25kGy

Complies with applicable parts of ISO 11137  
series

##### QUALITY ASSURANCE

Manufactured under a BS EN ISO  
13485:2016 quality management system

##### SHELF LIFE

3 years from date of manufacture

##### STORAGE CONDITIONS

Store in a dry, well ventilated area and  
protect from direct sunlight and heat  
sources

## TECHNICAL INFORMATION

Compatible with gamma irradiation  
BSE/TSE – considered not to present any risk\*  
Phthalate free  
Latex free  
PVC free  
Single scale easy to read graduation

\* raw materials may contain tallow-derived substances which have been produced under conditions which meet or exceed the requirements of EN 22442-1 and EMA 410/01 and are therefore considered not to present any risk in respect of TSE/BSE or other animal-borne diseases.

## ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
03ML05CC	3ml luer lock syringe plus tip protector, 5 pack	3ml	5	100	280 x 170 x 200mm	0.80 KG
03ML10CC	3ml luer lock syringe plus tip protector, 10 pack	3ml	10	100	280 x 170 x 200mm	0.55 KG
03ML25CC	3ml luer lock syringe plus tip protector, 25 pack	3ml	25	250	280 x 170 x 200mm	1.00 KG
03ML25CH	3ml luer lock syringe plus tip protector and cap, 25 pack	3ml	25	250	280 x 170 x 200 mm	1.00 KG
03ML25SC	3ml luer lock syringe plus barcode and cap, sleeved, 25 x 12 pack	3ml	25	300	280 x 340 x 400mm	3.30 KG
03ML50CC	3ml luer lock syringe plus tip protector, 50 pack	3ml	50	500	280 x 340 x 400mm	2.50 KG