



PRODUCT SPECIFICATION

STERILE PACKS

BD 3ML

SYRINGE MULTIPACKS

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Three piece BD single use syringe packs

PRODUCT FEATURES

Materials

Barrel and plunger – polypropylene
Gasket – latex-free, medical grade rubber
Lubricant – medical grade silicone oil
($<0.25 \text{ mg/cm}^2$)

Connection type

Luer-Lock (central)

Applicable standards

ISO 7886-1 / ISO 7886-2

Medical device classification

Class I (sterile)

MANUFACTURING

Becton Dickinson & Company 1 Becton
Drive Franklin Lakes NJ 07417 USA

PACKAGING

Four layer laminate
Easy tear notch
Alcohol and VHP compatible

STERILISATION DETAILS

Gamma irradiation min. 25kGy
Complies with applicable parts of ISO 11137
series

QUALITY ASSURANCE

Manufactured under a BS EN ISO
13485:2016 quality management system

SHELF LIFE

3 years from date of manufacture

STORAGE CONDITIONS

Store in a dry, well ventilated area and
protect from direct sunlight and heat
sources

TECHNICAL INFORMATION

Compatible with gamma irradiation
 BSE/TSE – considered not to present any risk*
 Phthalate free
 Latex free
 PVC free
 Single scale easy to read graduation

* raw materials may contain tallow-derived substances which have been produced under conditions which meet or exceed the requirements of EN 22442-1 and EMA 410/01 and are therefore considered not to present any risk in respect of TSE/BSE or other animal-borne diseases.

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
03ML05CC	3ml BD Syringe Sterile 20 Packs of 5	3ml	5	100	280 x 170 x 200mm	0.80 KG
03ML10CC	3ml BD Syringe Sterile 10 Packs of 10	3ml	10	100	280 x 170 x 200mm	0.55 KG
03ML25CC	3ml BD Syringe Sterile 10 Packs of 25	3ml	25	250	280 x 170 x 200mm	1.00 KG
03ML25CH	3ml BD Syringe and Cap Sterile 10 Packs of 25	3ml	25	250	280 x 170 x 200 mm	1.00 KG
03ML25SC	33ml BD Syringe and Cap Sleeved Sterile 12 Packs of 25	3ml	25	300	280 x 340 x 400mm	3.30 KG
03ML50CC	3ml BD Syringe Sterile 10 Packs of 50	3ml	50	500	280 x 340 x 400mm	2.50 KG