



PRODUCT SPECIFICATION

STERILE PACKS

BD 50ML

SYRINGE MULTIPACKS

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Three piece BD single use syringe packs

PRODUCT FEATURES

Materials

Barrel and plunger – polypropylene

Gasket – latex-free, medical grade rubber

Lubricant – medical grade silicone oil (<0.25 mg/cm²)

Connection type

Luer-Lock (central)

Applicable standards

ISO 7886-1

Medical device classification

Class I (sterile/measuring)

MANUFACTURING

Becton Dickinson & Company 1 Becton Drive
Franklin Lakes NJ 07417 USA

PACKAGING

Four layer laminate

Easy tear notch

Alcohol and VHP compatible

STERILISATION DETAILS

Gamma irradiation min. 25kGy

Complies with applicable parts of ISO 11137 series

QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016 quality management system

SHELF LIFE

3 years from date of manufacture

STORAGE CONDITIONS

Store in a dry, well ventilated area and protect from direct sunlight and heat sources

TECHNICAL INFORMATION

Compatible with gamma irradiation
 BSE/TSE – considered not to present any risk*
 Phthalate free
 Latex free
 PVC free
 Single scale easy to read graduation

* raw materials may contain tallow-derived substances which have been produced under conditions which meet or exceed the requirements of EN 22442-1 and EMA 410/01 and are therefore considered not to present any risk in respect of TSE/BSE or other animal-borne diseases.

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
50ML05CC	50ml BD Syringe Sterile 10 Packs of 5	50ml	5	50	280 x 340 x 400mm	2.30 KG
50ML10CC	50ml BD Syringe Sterile 10 Packs of 10	50ml	10	100	280 x 340 x 400mm	4.20 KG
50ML25CC	50ml BD Syringe Sterile 8 Packs of 25	50ml	25	200	560 x 330 x 380mm	8.00 KG
50ML25CH	50ml BD Syringe and Cap Sterile 8 Packs of 25	50ml	25	200	560 x 330 x 380mm	8.00 KG
50ML25SC	50ml BD Syringe and Cap Sleeved Sterile 6 Packs of 10	50ml	10	60	560 x 330 x 380mm	5.70 KG
50ML300C	50ml BD Syringe and Cap Sleeved Sterile 12 Packs of 25	50ml	25	300	550 x 500 x 420mm	14.20 KG