



PRODUCT SPECIFICATION

STERILE PACKS

BD 5ML

SYRINGE MULTIPACKS

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Three piece BD single use syringe packs

PRODUCT FEATURES

Materials

Barrel and plunger – polypropylene
Gasket – latex-free, medical grade rubber
Lubricant – medical grade silicone oil
($<0.25 \text{ mg/cm}^2$)

Connection type

Luer-Lock (central)

Applicable standards

ISO 7886-1 / ISO 7886-2
ISO 80369-7

Medical device classification

Class I (sterile)

MANUFACTURING

Becton Dickinson & Company 1 Becton
Drive Franklin Lakes NJ 07417 USA

PACKAGING

Four layer laminate
Easy tear notch
Alcohol and VHP compatible

STERILISATION DETAILS

Gamma irradiation min. 25kGy
Complies with applicable parts of ISO 11137
series

QUALITY ASSURANCE

Manufactured under a BS EN ISO
13485:2016 quality management system

SHELF LIFE

3 years from date of manufacture

STORAGE CONDITIONS

Store in a dry, well ventilated area and
protect from direct sunlight and heat
sources

TECHNICAL INFORMATION

Compatible with gamma irradiation
BSE/TSE – considered not to present any risk*
Phthalate free
Latex free
PVC free
Single scale easy to read graduation

*raw materials may contain tallow-derived substances which have been produced under conditions which meet or exceed the requirements of EN 22442-1 and EMA 410/01 and are therefore considered not to present any risk in respect of TSE/BSE or other animal-borne diseases.

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
05ML05CC	5ml luer lock syringe plus tip protector, 5 pack	5ml	5	100	280 x 170 x 200mm	0.90 KG
05ML10CC	5ml luer lock syringe plus tip protector, 10 pack	5ml	10	100	280 x 170 x 200mm	0.90 KG
05ML25CC	5ml luer lock syringe plus tip protector, 25 pack	5ml	25	250	280 x 170 x 390mm	1.70 KG
05ML25CH	5ml luer lock syringe plus tip protector and cap, 25 pack	5ml	25	250	280 x 170 x 390 mm	1.70 KG
05ML50CC	5ml luer lock syringe plus tip protector, 50 pack	5ml	50	500	280 x 340 x 400 mm	3.50 KG