



PRODUCT SPECIFICATION

DELTA TRIGGER SPRAY 950ml STERILE

CleanGuard **5**

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Neutral Detergent 950ml trigger spray.
Available in sterile version.
Safe to use on any surface and is colourless and odourless.
Pre-mixed solution ensures homogenous distribution.
Manufactured with WFI.
pH 6.5 - 7.5

PRODUCT FEATURES

Trigger spray with spray/off settings.
Trigger designed for low actuation force.
Tamper resistant fastening.
HDPE bottle.
A hydrophobic PTFE membrane preventing liquid and microbial penetration while allowing gas exchange in both directions to equalize pressure.

TECHNICAL INFORMATION

MANUFACTURING

The product is filled and capped in a Grade B (in use) cleanroom.
Packaging is performed in a Grade C (in use) cleanroom.
The neutral detergent is filtered through a 0.2 micron filter at the point of fill.

PACKAGING

Primary packaging: natural HDPE bottle with tamper-evident trigger spray.
Secondary packaging: double bagged in polythene pouches.

STERILISATION DETAILS

Gamma Irradiation (minimum 25kGy).

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001:2015 registered systems.

SHELF LIFE

2 years from date of manufacture.

STORAGE CONDITIONS

Store in a dry, well-ventilated area in ambient conditions.

ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
	ZSND05950	Delta 950ml Trigger Spray Sterile	950ml	Each	15 x 950ml	442 x 390 x 320mm	17 Kg

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please contact us.