



PRODUCT SPECIFICATION

DI Water 5 Litre RFU STERILE/NON-STERILE



PRODUCT INFORMATION

PRODUCT DESCRIPTION

Sterile and Non-sterile DI water in RFU format.

PRODUCT FEATURES

5 Litre container complete with a screw on tamper evident air tight cap.
Convenient carrying size and weight.

TECHNICAL INFORMATION

MANUFACTURING

The product is filled and capped in a Grade B (in use) cleanroom.

Packaging is performed in a Grade C (in use) cleanroom.

PACKAGING

Primary packaging: 5 Litre container complete with a screw on air tight cap.
Secondary packaging: double bagged in polythene bags.
2 double-bagged bottles in the box.

STERILISATION DETAILS

Gamma Irradiation (minimum 25kGy).

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems.

SHELF LIFE

2 years from date of manufacture.

STORAGE CONDITIONS

Store in a dry, well-ventilated area in ambient conditions.

ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
	ZSDI82005	DI Water 5L RFU Sterile	5L	Each	2 x 5L	335mm x 330mm x 265mm	11 Kg
ZNDI82005		DI Water 5L RFU Non-Sterile	5L	Each	2 x 5L	335mm x 330mm x 265mm	11 Kg