



## PRODUCT SPECIFICATION

### ECLIPSE SAFETY NEEDLES

#### PRODUCT INFORMATION

##### PRODUCT DESCRIPTION

Eclipse safety needles.  
20G x 40mm and 21G x 25mm.

##### PRODUCT FEATURES

###### Materials

Cannula - stainless steel.  
Hub, cover and shield - polypropylene.

###### Connection type

Luer lock.

###### Medical Device Classification

Not a medical device (pharmacy compounder).

##### MANUFACTURING

BD Spain.

##### PACKAGING

Four layer laminate.  
Easy tear notch.  
Alcohol and VHP compatible.  
Validated three year shelf-life.

##### STERILISATION DETAILS

Gamma irradiation min. 25kGy.  
Complies with applicable parts of  
ISO 11137 series.

##### QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016  
quality management system.

##### SHELF LIFE

3 years from date of manufacture.

##### STORAGE CONDITIONS

Store in a dry, well ventilated area and  
protect from direct sunlight and heat  
sources.

## TECHNICAL INFORMATION

Compatible with gamma irradiation.

BSE/TSE - considered not to present any risk.\*

Phthalate free.

Latex free.

\*Raw materials may contain tallow-derived substances which have been produced under conditions which meet or exceed the requirements of EN 22442-1 and EMA 410/01 and are therefore considered not to present any risk in respect of TSE/BSE or other animal-borne diseases.

## ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
20EC50CC-A	BD 20G Eclipse Needle Sterile 10 Packs of 50	n/a	50	500	325 x 280 x 390mm	2.40 KG
21EC50CC-A	BD 21G Eclipse Needle Sterile 10 Packs of 50	n/a	50	500	325 x 280 x 390mm	2.40 KG