



PRODUCT SPECIFICATION

ECLIPSE SAFETY NEEDLES

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Eclipse safety needles
20G x 40mm and 21G x 25mm

PRODUCT FEATURES

Materials

Cannula - stainless steel
Hub, cover and shield - polypropylene

Connection type

Luer lock

Medical Device Classification

Not a medical device (pharmacy compounder)

MANUFACTURING

BD Spain

PACKAGING

Four layer laminate
Easy tear notch
Alcohol and VHP compatible
Validated three year shelf-life

STERILISATION DETAILS

Gamma irradiation min. 25kGy
Complies with applicable parts of
ISO 11137 series

QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016
quality management system

SHELF LIFE

3 years from date of manufacture

STORAGE CONDITIONS

Store in a dry, well ventilated area and
protect from direct sunlight and heat
sources

TECHNICAL INFORMATION

Compatible with gamma irradiation
BSE/TSE free
Phthalate free
Latex free

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
20EC50CC-A	BD 20G Eclipse Needle Sterile 10 Packs of 50	n/a	50	500	325 x 280 x 390mm	2.40 KG
21EC50CC-A	BD 21G Eclipse Needle Sterile 10 Packs of 50	n/a	50	500	325 x 280 x 390mm	2.40 KG