



PRODUCT SPECIFICATION

FILTER NEEDLES STERILE PACKS

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Filter needles 19g x 1½".

PRODUCT FEATURES

Materials

Needle: stainless steel.
Cap: polyethylene.
Hub: polypropylene.

Connection type

Luer lock.

Medical device classification

Not a medical device (pharmacy compounder).

MANUFACTURING

B Braun, Malaysia.

PACKAGING

Four layer laminate.
Easy tear notch.
Alcohol and VHP compatible.
Validated three year shelf-life.

STERILISATION DETAILS

Gamma irradiation min. 25kGy.
Complies with applicable parts of ISO 11137 series.

QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016 quality management system.

SHELF LIFE

3 years from date of manufacture.

STORAGE CONDITIONS

Store in a dry, well ventilated area and protect from direct sunlight and heat sources.



TECHNICAL INFORMATION

Compatible with gamma irradiation.
BSE/TSE - considered not to present any risk.*
Phthalate free.
Latex free.

*This product contains materials of animal origin that are limited to tallow or tallow derivatives. The materials are rigorously processed according to prevailing guidance and standards in order to reduce or eliminate the risk of transmission of BSE/TSE.

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
19FN25CC	B Braun filter needle 19g x 1.5", 25 pack	n/a	25	250	290 x 170 x 200mm	0.80 kg
19FN25TB	B Braun filter needle 19g x 1.5", 25 pack, triple-wrapped	n/a	25	250	280 x 160 x 390mm	1.00 kg

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please contact us.