



PRODUCT SPECIFICATION

FILTER STRAW MULTIPACKS

PRODUCT INFORMATION

PRODUCT DESCRIPTION

B Braun filter straws

PRODUCT FEATURES

Materials

ABS / Nylon / DEHP-Free PVC

Connection type

Luer lock

Medical Device Classification

Not a medical device (pharmacy compounder)

MANUFACTURING

B Braun USA

PACKAGING

Four layer laminate

Easy tear notch

Alcohol and VHP compatible

Validated three year shelf-life

STERILISATION DETAILS

Gamma irradiation min. 25kGy

Complies with applicable parts of ISO 11137 series

QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016 quality management system

SHELF LIFE

3 years from date of manufacture

STORAGE CONDITIONS

Store in a dry, well ventilated area and protect from direct sunlight and heat sources

TECHNICAL INFORMATION

Compatible with gamma irradiation
BSE/TSE free
Phthalate free
Latex free

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
FRSW25CC	Filter Straw Sterile 10 Packs of 25	n/a	25	250	290 x 170 x 200mm	0.66 KG