



PRODUCT SPECIFICATION

BLACK CLEANROOM MARKER PEN STERILE

GuardMark 

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Sterile individually packed marker pens.

PRODUCT FEATURES

New improved Micronclean branded marker pen with quick-drying ink which is, water and UV light resistant. Although the ink is quick drying our new pen does not dry out even if the cap is left off for a day.

Our new pen has a proven ability to write on most typical laboratory surfaces including plastic, glass, and stainless steel.

In addition, the ink in our new pen does not dry out even if the cap is left off for a day.

TECHNICAL INFORMATION

PACKAGING

Primarily packaged in polythene bags and secondary packaged in a cardboard carton in a Grade C (in use) Cleanroom.

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems.

STERILISATION DETAILS

Gamma Irradiation (minimum 25kGy).

SHELF LIFE

2 year from date of manufacture.

STORAGE CONDITIONS

Store in ambient conditions.

ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
	ZSMK10110	GuardMark 2 Cleanroom Marker Pen Black Sterile	134x17mm	1	10	190 x 101 x 51mm	0.16 Kg