



PRODUCT SPECIFICATION

IPA MICROFIBRE MOP STERILE

PureGuard **1**

PRODUCT INFORMATION

PRODUCT DESCRIPTION

White single use microfibre sterile disposable mop pre-saturated with 70% IPA / 30%WFI designed for hard surface mopping within a cleanroom environment.

Uniformly dosed to achieve 15m² coverage eliminating the need for in cleanroom dosing or mixing.

PRODUCT FEATURES

Materials

Mop front made from Microfibre (80% polyester + 20% polyamide)

Mop back made from 100% polyester
2 pocket design enables hands free and secure connection to the mop frame
Mop size 440mm x 140mm

MANUFACTURING

Washed within a ISO 14644-1 Class 4 (in use) Cleanroom laundry

The product is pre-saturated within an ISO 14644-1 Class 5 (in use) unidirectional airflow cabinets

Packaging is performed within an ISO 14644-1 Class 7 (in use) cleanroom

IPA is filtered through a 0.2 micron filter at the point of pre-saturation

PACKAGING

Primary packaging: double bagged in polythene bags

Secondary packaging: cardboard box

STERILISATION DETAILS

Gamma Irradiation (minimum 25kGy)

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems.

SHELF LIFE

2 years from date of manufacture

STORAGE CONDITIONS

Store in a dry, well-ventilated area in ambient conditions.



ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
MSMI21016 - A	PureGuard 1 IPA Multi Packed Disposable Microfibre Mop Sterile	440mm x 140mm	10	3 x 10	450 x 200 x 200mm	8.5 Kg
MSMI24016 - A	PureGuard 1 IPA Individually Packed Disposable Microfibre Mop Sterile	440mm x 140mm	Each	30 x 1	450 x 200 x 200mm	8 Kg