



## PRODUCT SPECIFICATION

### BETA MICROFIBRE MOP STERILE

PureGuard **3**

#### PRODUCT INFORMATION

##### PRODUCT DESCRIPTION

White single use microfibre sterile disposable mop pre-saturated with Beta disinfectant designed for hard surface mopping within a cleanroom environment.

Uniformly dosed to achieve 15m<sup>2</sup> coverage eliminating the need for in cleanroom dosing or mixing.

##### PRODUCT FEATURES

###### Materials

Mop front made from Microfibre (80% polyester + 20% polyamide)

Mop back made from 100% polyester

2 pocket design enables hands free and secure connection to the mop frame

Mop size 440mm x 140mm

##### MANUFACTURING

Washed within a ISO 14644-1 Class 4 (in use) Cleanroom laundry

The product is pre-saturated within an ISO 14644-1 Class 5 (in use) unidirectional airflow cabinets

Packaging is performed within an ISO 14644-1 Class 7 (in use) cleanroom

Beta is filtered through a 0.2 micron filter at the point of pre-saturation

##### PACKAGING

Primary packaging: double bagged in polythene bags

Secondary packaging: cardboard box

##### STERILISATION DETAILS

Gamma Irradiation (minimum 25kGy)

##### QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems

##### SHELF LIFE

2 years from date of manufacture

##### STORAGE CONDITIONS

Store in a dry, well-ventilated area in ambient conditions



## ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
MSMB24016-A	PureGuard 3 Beta Microfibre Mop Multi Packed Sterile	440mm x 140mm	10	3 x 10	450 x 200 x 200mm	8.5 Kg
MSMB21016-A	PureGuard 3 Beta Microfibre Mop Individually Packed Sterile	440mm x 140mm	Each	30 x 1	450 x 200 x 200mm	8 Kg