



PRODUCT SPECIFICATION

ALPHA MICROFIBRE MOP STERILE

PureGuard **4**

PRODUCT INFORMATION

PRODUCT DESCRIPTION

White single use microfiber sterile disposable mop pre-saturated with Alpha disinfectant designed for hard surface mopping within a cleanroom environment.

Uniformly dosed to achieve 15m² coverage eliminating the need for in cleanroom dosing or mixing.

PRODUCT FEATURES

Materials

Mop front made from Microfibre (80% polyester + 20% polyamide).

Mop back made from 100% polyester.

2 pocket design enables hands free and secure connection to the mop frame.

Mop size 440mm x 140mm.

TECHNICAL INFORMATION

MANUFACTURING

The product is filled and sealed in a Grade B (in use) cleanroom.

Packaging is performed in a Grade C (in use) cleanroom.

The disinfectant is filtered through a 0.2 micron filter at the point of fill.

PACKAGING

Primary packaging: double bagged in polythene bags.

Secondary packaging: cardboard box.

STERILISATION DETAILS

Gamma Irradiation (minimum 25kGy).

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems.

SHELF LIFE

2 years from date of manufacture.

STORAGE CONDITIONS

Store in a dry, well-ventilated area in ambient conditions.



ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
MSMA24016	PureGuard 4 Alpha Microfibre Mop Individually Packed Sterile	440mm x 140mm	Each	30 x 1	450 x 200 x 200mm	8 Kg

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please contact us.