



## PRODUCT SPECIFICATION

### DELTA MICROFIBRE MOP STERILE

PureGuard **5**

#### PRODUCT INFORMATION

##### PRODUCT DESCRIPTION

White single use microfibre sterile disposable mop pre-saturated with Delta neutral detergent designed for hard surface mopping within a cleanroom environment.

Uniformly dosed to achieve 15m<sup>2</sup> coverage eliminating the need for in cleanroom dosing or mixing.

##### PRODUCT FEATURES

###### Materials

Mop front made from Microfibre (80% polyester + 20% polyamide)

Mop back made from 100% polyester  
2 pocket design enables hands free and secure connection to the mop frame

Mop size 440mm x 140mm

##### MANUFACTURING

Washed within a ISO 14644-1 Class 4 (in use) Cleanroom laundry.

The product is pre-saturated within an ISO 14644-1 Class 5 (in use) unidirectional airflow cabinets.

Packaging is performed within an ISO 14644-1 Class 7 (in use) cleanroom

Neutral Detergent is filtered through a 0.2 micron filter at the point of pre-saturation.

##### PACKAGING

Primary packaging: double bagged in polythene bags

Secondary packaging: cardboard box

##### STERILISATION DETAILS

Gamma Irradiation (minimum 25kGy)

##### QUALITY ASSURANCE

Product supplied under BS EN ISO 9001:2015 registered systems.

##### SHELF LIFE

2 years from date of manufacture

##### STORAGE CONDITIONS

Store in a dry, well-ventilated area in ambient conditions.

##### HAZARD STATEMENT

Non-hazardous for full details please consult the relevant SDS.



## ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
MSMN21016-A	PureGuard 5 Delta Multi Packed Disposable Microfibre Mop Sterile	440mm x 140mm	10	3 x 10	450 x 200 x 200mm	8.5 Kg
MSMN24016-A	PureGuard 5 Delta Individually Packed Disposable Microfibre Mop Sterile	440mm x 140mm	Each	30	450 x 200 x 200mm	8.5 Kg