



PRODUCT SPECIFICATION

DUAL LAYER POLYESTER ISOLATOR BONNET STERILE

PureGuard **6**

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Universal fit isolator cleaning tool made from dual layer continuous filament polyester wipes designed specifically to help clean inside of isolators and LAF cabinets.

PRODUCT FEATURES

Low linting, dual layer 100% continuous filament polyester for increased absorbency;
Universal fit;
Supplied sterile one per bag ready for use;
Aqueous processing at ISO 14644-1 EU GMP Grade B.

TECHNICAL INFORMATION

MANUFACTURING

Aqueous processing at GMP Grade B.

PACKAGING

Primarily packaged in polythene bags and secondary packaged in a cardboard carton in a Grade C (in use) Cleanroom.

STERILISATION DETAILS

Gamma Irradiation (minimum 25kGy).

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems.

SHELF LIFE

2 years from date of manufacture;
Use on day of opening.

STORAGE CONDITIONS

Store in cool, well ventilated area.

PRODUCT SPECIFICATION

Material	100% continuous filament polyester
Weight	250g/m ²
Particles by sizes @ ≥ 0.5um (IEST-RP-CC003.23), Helmke drum @ 10rpm	≤ 300 counts Class 1

ORDERING/PACKING INFORMATION

Sterile Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
ISED0140	Dry Dual Layer Polyester Isolator Bonnet Sterile	One Size	Each	40	245 x 195 x 260mm	1 Kg