



# PRODUCT SPECIFICATION

## DUAL LAYER POLYESTER ISOLATOR BONNET STERILE

PureGuard **6**

### PRODUCT INFORMATION

#### PRODUCT DESCRIPTION

Universal fit isolator cleaning tool made from dual layer continuous filament polyester wipes designed specifically to help clean inside of isolators and LAF cabinets.

#### PRODUCT FEATURES

Low linting, dual layer 100% continuous filament polyester for increased absorbency;  
Universal fit;  
Supplied sterile one per bag ready for use;  
Aqueous processing at ISO 14644-1 EU GMP Grade B;  
Tested to ASTM-D-629.

#### MANUFACTURING

Aqueous processing at ISO 14644-1 EU GMP Grade B;  
Packaging is performed within an ISO 14644-1 Class 7 (in use) cleanroom.

#### PACKAGING

Primary packaging: Double bagged in polythene packaging;  
Supplied 40 individually packed bonnets per carton.

#### STERILISATION DETAILS

Gamma Irradiation (minimum 25kGy)

#### QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems.

#### SHELF LIFE

2 years from date of manufacture;  
Use on day of opening.

#### STORAGE CONDITIONS

Store in cool, well ventilated area.

#### MATERIAL DIMENSIONS

Wipe Material: 100% continuous filament polyester.  
Wipe Weight: 250g/m<sup>2</sup>

### ORDERING/PACKING INFORMATION

Sterile Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
ISED0140	Dry Dual Layer Polyester Isolator Bonnet Sterile	One Size	Each	40	245 x 195 x 260mm	1 Kg

ISSUE 1 | DATE 12/03/2021 | PAGE 1 OF 1