



PRODUCT SPECIFICATION

STERICAN NEEDLES STERILE PACKS

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Compounding needles 19g & 21g x 1½"

PRODUCT FEATURES

Materials

Needle: stainless steel
Cap: polyethylene
Hub: polypropylene

Connection type

Luer lock

Medical device classification

Not a medical device (pharmacy compounder)

MANUFACTURING

B Braun, Malaysia

PACKAGING

Four layer laminate
Easy tear notch
Alcohol and VHP compatible
Validated three year shelf-life

STERILISATION DETAILS

Gamma irradiation min. 25kGy
Complies with applicable parts of
ISO 11137 series

QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016
quality management system

SHELF LIFE

3 years from date of manufacture

STORAGE CONDITIONS

Store in a dry, well ventilated area and
protect from direct sunlight and heat
sources

TECHNICAL INFORMATION

Compatible with gamma irradiation
BSE/TSE free
Phthalate free
Latex free



ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
19GN50CC	B Braun needle 19g x 1.5", 50 pack	n/a	50	500	290 x 170 x 200mm	0.95 kg
21GN50CC	B Braun needle 21g x 1.5", 50 pack	n/a	50	500	290 x 170 x 200mm	0.95 kg
19GN50TB	B Braun needle 19g x 1.5", 50 pack, triple-wrapped	n/a	50	500	290 x 170 x 200mm	0.95 kg
21GN50TB	B Braun needle 21g x 1.5", 50 pack, triple-wrapped	n/a	50	500	290 x 170 x 200mm	0.95 kg

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please contact us.