



PRODUCT SPECIFICATION

TERUMO NEEDLES STERILE PACKS

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Needles for pharmacy compounding
18g, 19g & 21g x 40mm

PRODUCT FEATURES

Materials

Needle: stainless steel
Cap: polypropylene
Hub: polypropylene
Sealing label: self-adhesive paper (tamper evident closure)

Connection type

Luer lock

Medical device classification

Not a medical device (pharmacy compounder)

MANUFACTURING

Micronclean are the legal manufacturer of this product and take full responsibility for its safety and performance

This product is intended for pharmacy compounding use only

The CE mark 0197 and reference to Terumo Europe as legal manufacturer present on the individual needles of this product should be disregarded

PACKAGING

Four layer laminate
Easy tear notch
Alcohol and VHP compatible
Validated three year shelf-life

STERILISATION DETAILS

Gamma irradiation min. 25kGy
Complies with applicable parts of ISO 11137 series

QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016 quality management system

SHELF LIFE

3 years from date of manufacture

STORAGE CONDITIONS

Store in a dry, well ventilated area and protect from direct sunlight and heat sources

TECHNICAL INFORMATION

Compatible with gamma irradiation
BSE/TSE free
Phthalate free
Latex free

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
18TR50CC	Terumo Needle 18g x 40mm Sterile	n/a	50	500	290 x 170 x 200mm	1.00 KG
19TR50CC	Terumo Needle 19g x 40mm Sterile	n/a	50	500	290 x 170 x 200mm	1.00 KG
21TR50CC	Terumo Needle 21g x 40mm Sterile	n/a	50	500	290 x 170 x 200mm	0.80 KG