



PRODUCT SPECIFICATION

BLUE POLYTHENE OVERSHOES NON-STERILE

TuffGuard **2**

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Blue disposable heavyweight overshoe

PRODUCT FEATURES

Elastic opening for easy donning
Lightweight
Durable
Cost effective

PACKAGING

Primary packaged in polythene packaging

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems

STORAGE CONDITIONS

Store in a dry, well-ventilated area in ambient conditions

MATERIAL AND PERFORMANCE

LDPE 35%; LLDPE 65% (polyethylene) material
Material = 40µm embossed
14" Overshoe size = 360mm ±5% length
150mm ±5% width
16" Overshoe size = 410mm ±5% length
150mm ±5% width

ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Case Qty
SNBL31014		Blue Polythene Overshoe 14" Non-Sterile	360mm x 150mm	100	20 x 100
SNBL31016		Blue Polythene Overshoe 16" Non-Sterile	410mm x 150mm	100	20 x 100