



PRODUCT SPECIFICATION

ALPHA SESSIONAL POLYCELLULOSE C-FOLDED POUCH WIPES STERILE



PRODUCT INFORMATION

PRODUCT DESCRIPTION

Alpha Pre-saturated non-woven Polycellulose sterile wipes
Supplied C-folded in a pouch ready immediate use
Increased disinfectant and retention from re-sealable pouch
Specifically designed for use in Cleanrooms

PRODUCT FEATURES

Non-woven polycellulose wipe with sealed edges
Pre-saturated with Alpha to ensure optimum disinfection
Supplied in a small pack size to minimise waste

TECHNICAL INFORMATION

MANUFACTURING

The product is filled and sealed in a Grade B (in use) cleanroom. Packaging is performed in a Grade C (in use) cleanroom.
The disinfectant is filtered through a 0.2 micron filter at the point of fill.

PACKAGING

Supplied double bagged.
Secondary Packaging: Cardboard carton.

STERILISATION DETAILS

Sterilised by Gamma Irradiation (minimum 25kGy)

SHELF LIFE

Two years from date of manufacture

STORAGE CONDITIONS

Store in a cool dry place (< 40°C) away from direct sunlight and fluorescent light

PRODUCT INFORMATION

DIMENSIONS

Wipe Size: 23cm x 23cm (9" x 9")

DISINFECTANT PERFORMANCE

SUSPENSION TESTS							
	Protocol		Organism	Contact Time (min)	Conditions	Log reduction to pass	Log reduction achieved
Suspension Tests Phase 2 Step 1	EN 1276	Bactericidal	Pseudomonas aeruginosa	5 mins	Dirty	> 5.0	7.76
			Escherichia coli	5 mins	Dirty	> 5.0	7.40
			Staphylococcus aureus	5 mins	Dirty	> 5.0	7.72
			Enterococcus hirae	5 mins	Dirty	> 5.0	7.83
	EN 1650	Fungicidal/ Yeasticidal	Aspergillus niger	30 mins	Dirty	> 4.0	> 4.00
			Candida albicans	30 mins	Dirty	> 4.0	> 4.00
Suspension Tests Phase 2 Step 2	EN 13697	Bactericidal	Pseudomonas aeruginosa	5 mins	Clean	> 4.0	4.42
			Escherichia coli	5 mins	Clean	> 4.0	6.08
			Staphylococcus aureus	5 mins	Clean	> 4.0	5.04
			Enterococcus hirae	5 mins	Clean	> 4.0	4.39
			Corynebacterium diptheriae	5 mins	Clean	> 4.0	4.99

ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Pack Qty	Case Qty
	WSAP01009	Alpha Sessional Polycellulose C-folded Pouch Wipes Sterile	230 x 230mm	15 wipes per Pouch	15	10 x 15

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please contact us.