



PRODUCT SPECIFICATION

DELTA SESSIONAL POLYCELLULOSE C-FOLDED POUCH WIPES STERILE



PRODUCT INFORMATION

PRODUCT DESCRIPTION

Delta Pre-saturated non-woven Poly-cellulose sterile wipes.
Supplied C-folded in a pouch ready immediate use.
Increased disinfectant and retention from re-sealable pouch.
Specifically designed for use in Cleanrooms.

PRODUCT FEATURES

Non-woven polycellulose wipe with sealed edges.
Pre-saturated with Delta to ensure optimum cleaning.
Supplied in a small pack size to minimise waste.

TECHNICAL INFORMATION

MANUFACTURING

The product is filled and sealed in a Grade B (in use) cleanroom. Packaging is performed in a Grade C (in use) cleanroom.
The disinfectant is filtered through a 0.2 micron filter at the point of fill.

PACKAGING

Supplied double bagged.
Secondary Packaging: Cardboard carton.

STERILISATION DETAILS

Sterilised by Gamma Irradiation (minimum 25kGy).

SHELF LIFE

Two years from date of manufacture.

STORAGE CONDITIONS

Store in a cool dry place (< 40°C) away from direct sunlight and fluorescent light.

PRODUCT INFORMATION

DIMENSIONS

Wipe Size: 23cm x 23cm (9" x 9")

ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Pack Qty	Case Qty
	WSND05009	Delta Sessional Poly-cellulose C-folded Pouch Wipes Sterile	230 x 230mm	15 wipes per Pouch	15	10 x 15

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please contact us.