



PRODUCT SPECIFICATION

DELTA SESSIONAL POLY-CELLULOSE C-FOLDED POUCH WIPES STERILE



PRODUCT INFORMATION

PRODUCT DESCRIPTION

Delta Pre-saturated non-woven Poly-cellulose sterile wipes
Supplied C-folded in a pouch ready immediate use
Increased disinfectant and retention from re-sealable pouch
Specifically designed for use in Cleanrooms

PRODUCT FEATURES

Non-woven polycellulose wipe with sealed edges
Pre-saturated with Delta to ensure optimum cleaning
Supplied in a small pack size to minimise waste

MANUFACTURING

The product is pre-saturated in ISO 14644-1 Class 5 (in use) unidirectional airflow cabinets
The Neutral detergent is filtered through a 0.2 micron filter at the point of fill

PACKAGING

Supplied double bagged
All packaging is completed within an ISO 14644-1 Class 7 (in use) cleanroom

STERILISATION DETAILS

Sterilised by Gamma Irradiation (minimum 25kGy)

SHELF LIFE

Two years from date of manufacture

STORAGE CONDITIONS

Store in a cool dry place (< 40°C) away from direct sunlight and fluorescent light

PRODUCT INFORMATION

DIMENSIONS

Wipe Size: 23cm x 23cm (9" x 9")

ORDERING/PACKING INFORMATION

| Non-Sterile Code | Sterile Code | Description | Size | Pack Qty | Pack Qty | Case Qty |
|------------------|---------------|---|-------------|--------------------|----------|----------|
| | WSND05009-A/I | Delta Sessional Poly-cellulose C-folded Pouch Wipes Sterile | 230 x 230mm | 15 wipes per Pouch | 15 | 10 x 15 |