



# PRODUCT SPECIFICATION

## POLYESTER POUCH WIPE STERILE AND NON-STERILE



### PRODUCT INFORMATION

#### PRODUCT DESCRIPTION

VeriGuard 3 is a pre-saturated Polyester wipe designed for use in high grade Cleanrooms as a hard surface wipe.

#### PRODUCT FEATURES

100% continuous filament pouch wipe pre-saturated with IPA 70/ 30 WFI supplied in a resealable pouch.

Each wipe is 300mm x 300mm supplied either non-folded or V-Folded to suit your requirements.

#### MANUFACTURING

Primary packaged in polythene bags in ISO 14644 Class 5 (at rest) cleanroom.

#### DISINFECTANT SPECIFICATION

Contains: 70% Propan-2-ol (#67-63-0) and 30% Water for Injection grade water.

#### PACKAGING

Packed in cleanroom compatible pouch containing either 50 or 15 wipes.  
Double Bagged (sterile only).

#### QUALITY ASSURANCE

Product supplied under BS EN ISO 9001.

#### SHELF LIFE

2 year from date of manufacture.

#### STORAGE CONDITIONS

Store in a cool, dry, well ventilated area and protect from direct sunlight heat and naked flames.

## MATERIALS & PERFORMANCE

Performance Data			
Product Code	Weight	Size	Thickness
MCZNLA50N-A/I	130±10g/m <sup>2</sup>	300mm x 300mm	
MCZSLA50S-A/I	130±10g/m <sup>2</sup>	300mm x 300mm	
WNVG06030-A/I	120 + 5 gsm	300mm x 300mm	0.52±0.04mm
WSVG06030-A/I	120 + 5 gsm	300mm x 300mm	0.52±0.04mm

## PRODUCT INFORMATION

### DISINFECTANT PERFORMANCE

SUSPENSION TESTS					
Protocol	Organism	Contact Time (min)	Conditions	Log reduction to pass	Log reduction
BS EN 1040	S. aureus	5 min	Dirty	Log 5.0	> 5.3
	P. aeruginosa	5 min	Dirty	Log 5.0	> 5.4
BS EN 1275	C. albicans	5 min	Dirty	Log 4.0	> 4.4
BS EN 1276	S. aureus	5 min	Clean	Log 5.0	> 5.3
	E. coli	5 min	Clean	Log 5.0	> 5.3
	P. aeruginosa	5 min	Clean	Log 5.0	> 5.4
	E. hirae	5 min	Clean	Log 5.0	> 5.3
BS EN 1650	C. albicans	5 min	Clean	Log 4.0	> 4.4

SURFACE TESTS					
Protocol	Organism	Contact Time (min)	Conditions	Log reduction to pass	Log reduction
BS EN 13697	S. aureus	5 min	Clean	Log 4.0	6.0
	E. coli	5 min	Clean	Log 4.0	> 4.8
	P. aeruginosa	5 min	Clean	Log 4.0	> 4.1
	E. hirae	5 min	Clean	Log 4.0	> 5.5
	C. albicans	5 min	Clean	Log 3.0	> 4.0
	A. niger	5 min	Clean	Log 3.0	3.7

Please note: Data shown above are just typical values.

## ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
MCZNLA50N-A/I		IPA Polyester Non-folded Pouch Wipe Non-sterile	300mm x 300mm	50	5 x 50	390 x 390 x 190mm	8.0kg
	MCZSLA50S-A/I	IPA Polyester Non-folded Pouch Wipe Sterile	300mm x 300mm	50	5 x 50	390 x 390 x 190mm	8.0kg
WNVG06030-A/I		IPA Sessional Polyester V-folded Pouch Wipe Non-sterile	300mm x 300mm	15	30 x 15	400 x 400 x 395mm	15.78kg
	WSVG06030-A/I	IPA Sessional Polyester V-folded Pouch Wipe Sterile	300mm x 300mm	15	30 x 15	400 x 400 x 395mm	15.78kg