



PRODUCT SPECIFICATION

BETA SESSIONAL POLYCELLULOSE C-FOLDED POUCH WIPES STERILE



PRODUCT INFORMATION

PRODUCT DESCRIPTION

Beta Pre-saturated non-woven Polycellulose sterile wipes
Supplied C-folded in a pouch ready immediate use
Increased disinfectant and retention from re-sealable pouch
Specifically designed for use in Cleanrooms

PRODUCT FEATURES

Non-woven polycellulose wipe with sealed edges
Pre-saturated with Beta to ensure optimum disinfection
Supplied in a small pack size to minimise waste

MANUFACTURING

The product is pre-saturated in ISO 14644-1 Class 5 (in use) unidirectional airflow cabinets
The disinfectant is filtered through a 0.2 micron filter at the point of fill

PACKAGING

Supplied double bagged
All packaging is completed within an ISO 14644-1 Class 7 (in use) cleanroom

STERILISATION DETAILS

Sterilised by Gamma Irradiation (minimum 25kGy)

SHELF LIFE

Two years from date of manufacture

STORAGE CONDITIONS

Store in a cool dry place (< 40°C) away from direct sunlight and fluorescent light

PRODUCT INFORMATION

DIMENSIONS

Wipe Size: 23cm x 23cm (9" x 9")

DISINFECTANT PERFORMANCE

SUSPENSION TESTS							
	Protocol		Organism	Contact Time (min)	Conditions	Log reduction to pass	Log reduction achieved
Suspension Tests Phase 2 Step 1	EN 1276	Bactericidal	Salmonella typhimurium	1 min	Dirty	> 5.0	6.30
			Escherichia coli	1 min	Dirty	> 5.0	6.20
			Enterococcus faecalis	1 min	Dirty	> 5.0	7.30
			Enterococcus hirae	1 min	Dirty	> 5.0	6.60
	EN 1650	Fungicidal/ Yeasticidal	Aspergillus niger	2 mins	Clean	> 4.0	> 4.00
			Candida albicans	2 mins	Clean	> 4.0	> 4.00
	EN 13704	Sporicidal	Clostridium difficile spores	2 mins	Clean	> 3.0	> 3.00
Suspension Tests Phase 2 Step 2	EN 13697	Bactericidal	Pseudomonas aeruginosa	5 mins	Clean	> 4.0	5.70
			Escherichia coli	5 mins	Clean	> 4.0	> 6.56
			Staphylococcus aureus	5 mins	Clean	> 4.0	6.02
			Enterococcus hirae	5 mins	Clean	> 4.0	4.38
		Fungicidal	Aspergillus niger	15 mins	Dirty	> 3.0	4.08

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Pack Qty	Case Qty
WSBT02009 - A/I	Beta Sessional Polycellulose C-folded Pouch Wipes Sterile	230 x 230mm	15 wipes per Pouch	15	10 x 15