



PRODUCT SPECIFICATION

CHEMO DISPENSING PINS

STERILE PACKS

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Vented chemo dispensing pin for injecting or withdrawing cytostatic drugs from vials or semi-rigid containers.

Designed to eliminate aerosoling of chemotherapeutic medications during the reconstitution and dispensing process

PRODUCT FEATURES

Materials

Polyethylene / ABS

Connection type

Luer lock

Medical device classification

Not a medical device (pharmacy compounder)

MANUFACTURING

B Braun Medical Inc, USA

PACKAGING

Four layer laminate

Easy tear notch

Alcohol and VHP compatible

Validated three year shelf-life

STERILISATION DETAILS

Gamma irradiation min. 25kGy

Complies with applicable parts of

ISO 11137 series

QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016 quality management system

SHELF LIFE

3 years from date of manufacture

STORAGE CONDITIONS

Store in a dry, well ventilated area and protect from direct sunlight and heat sources



TECHNICAL INFORMATION

Compatible with gamma irradiation
BSE/TSE free
Phthalate free
Latex free

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
CDPB25CC	B Braun chemo dispensing pins, 25 pack	n/a	25	100	330 x 280 x 400mm	1,20 KG