



PRODUCT SPECIFICATION

DELTA PLUS 5 LITRE RFU STERILE



PRODUCT INFORMATION

PRODUCT DESCRIPTION

Delta Plus is a sterile, low-foaming, neutral detergent formulated for use in cleanrooms and other controlled environments.

Manufactured exclusively from naturally derived ingredients and blended with Water for Injection (WFI), it contains no petrochemical components and carries no hazard pictograms, signal words or hazard statements.

Designed to meet the requirements of EU GMP Annex 1 Section 4.33, Delta Plus provides effective removal of surface contamination, disinfectant residues and bioburden, offering a sustainable and user-friendly solution for the first stage of a validated two-step cleaning and disinfection programmes.

PRODUCT FEATURES

Hard surface cleaning agent.
5 Litre container complete with a screw on air tight cap.
Convenient carrying size and weight.
Visible liquid level.
Easy to use.

TECHNICAL INFORMATION

MANUFACTURING

The product is filled and capped in a Grade B (in use) cleanroom.

Packaging is performed in a Grade C (in use) cleanroom.

The disinfectant is filtered through a 0.2 micron filter at the point of fill.

PACKAGING

Primary packaging: 5 Litre container complete with a screw on air tight cap.
Secondary packaging: double bagged in polythene bags.
2 double-bagged bottles in the box.

STERILISATION DETAILS

Gamma Irradiation (minimum 25kGy).

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001:2015 registered systems.

SHELF LIFE

2 years from date of manufacture.

STORAGE CONDITIONS

Store in a dry, well-ventilated area in ambient conditions.

RECYCLING INFORMATION

Component Name:	Component Material:	Disposal:	Guidance:
5L Container	HDPE	Recyclable	Ensure Cap is Removed
Venting Bottle Cap	Polypropylene	Recyclable	Check With Local Authority, Remove From Container
Polybag	LDPE	Recyclable	Check With Local Authority, Remove Any Labels
Box	Cardboard	Recyclable	Remove Irradiation Labels
Product Label	Gloss Clear Laminated Polypropylene Label with an Acrylic Water-Borne Adhesive	Non-recyclable	Remove From Container
Box Label	Paper Based with Acrylic Water-Borne Adhesive	Non-recyclable Individually	Leave On Cardboard Box
Batch Label (Vellum Label)	Uncoated Paper with an Acrylic Adhesive	Non-recyclable Individually	Leave On Cardboard Box
Cardboard Reinforced Tape	Cardboard	Recyclable	Leave On Cardboard Box

ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
	RSND5LT02	Delta Plus 5L RFU Sterile	5L	Each	2 x 5L	290 x 215 x 290mm	10.5 KG

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please contact us.