



PRODUCT SPECIFICATION

CODAN DISPENSING PINS STERILE PACKS

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Dispensing pin for safe injection or withdrawal of cytostatic drugs from vials

PRODUCT FEATURES

Materials

ABS (housing)
Polyethylene (vent cap and sheath)

Connection type

Luer lock

Medical device classification

Not a medical device (pharmacy compounder)

MANUFACTURING

Codan, Germany

PACKAGING

Four layer laminate
Easy tear notch
Alcohol and VHP compatible
Validated three year shelf-life

STERILISATION DETAILS

Gamma irradiation min. 25kGy
Complies with applicable parts of ISO 11137 series

QUALITY ASSURANCE

Manufactured under a BS EN ISO 13485:2016 quality management system

SHELF LIFE

3 years from date of manufacture

STORAGE CONDITIONS

Store in a dry, well ventilated area and protect from direct sunlight and heat sources



TECHNICAL INFORMATION

Compatible with gamma irradiation
BSE/TSE free
Phthalate free
Latex free

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
COSP25CC	Codan dispensing pins, 25 pack	n/a	25	250	290 x 170 x 200mm	1.20 kg
COSP25TB	Codan dispensing pins, 25 pack, triple-wrapped	n/a	25	250	290 x 170 x 200mm	1.20 kg

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please contact us.